

# Allergy diagnostics

# Allergy diagnostics



# Specific IgE

Inhalation · Food · Atopy · Insect venoms



For more information on this subject scan the QR code or enter the Quick Link code q002 at [www.euroimmun.com](http://www.euroimmun.com)

## Inhalation

**Clinical information:** An allergy is a misguided reaction of the immune system to foreign substances that should normally be tolerated. Alongside any genetic predisposition, numerous non-genetic factors also play a role, such as exposure to the allergen, nutritional condition, existing chronic diseases and acute viral infections. The most frequently occurring allergy is a type I hypersensitivity reaction, in which specific IgE antibodies are formed. The symptoms (rash, oedema or itching) generally occur shortly after contact with the allergen. These allergies are therefore also termed immediate type reactions. More than 15% of the population of industrially developed countries suffer from allergic reactions of the immediate type. Typical allergic reactions are rhinitis, conjunctivitis and allergic asthma. An increase in allergic rhinitis has been observed worldwide with a prevalence of 10 to 40% in certain regions.

In the case of inhalation allergies, the allergens enter the body through the air and through the mucous membranes. Seasonal allergens (pollen from trees, grasses and herbs) play a role, as do also indoor allergens (house dust mites, domestic animals and mould spores) which occur the whole year round. The allergic symptoms intensify with every further exposure to the allergen. If a systemic allergic reaction occurs, serious, even life-threatening reactions can result (anaphylactic shock).



**Diagnostics:** A confirmed allergy diagnosis demands an extensive patient anamnesis and various diagnostic methods (in vivo and in vitro tests), depending on the patient's clinical symptoms. In individual cases, provocation tests can be additionally performed to help with the diagnosis. The results of the various diagnostic methods are of equal value and complement each other. All individual results must always be interpreted within the context of the anamnesis and the clinical findings.

Various **inhalation profiles** are available for the clarification of inhalation allergies. These permit the semi-quantitative and monospecific in vitro determination of human IgE antibodies against the most frequent inhalation allergens in serum or in plasma. Moreover, **country-specific inhalation profiles** are available which vary in their allergen composition and are optimised with regard to the regional conditions.

In addition, the total IgE concentration in the serum can be determined using the Total IgE ELISA to allow differentiation between allergic and intrinsic asthma, between allergic and vasomotor rhinitis and between atopic and seborrhoeic dermatitis.

## Product overview

Method	Substrate	Application	Order number	Page
EUROLINE	Inhalation (g1, g3, g6, g12, t2, t3, t4, t7, w1, w6, w9, d1, d2, e1, e2, e3, m1, m2, m3, m6, CCD)	Profile with the most important inhalation allergens, covering around 95% of all inhalative sensitisations	DP 3110-1601 E	209
	Paediatric Inhalation (g6, g12, t2, t3, t4, w6, w8, w9, d1, d2, e1, e2, e3, e6, e82, e84, m1, m2, m3, m6, CCD)	Extended inhalation profile with the most important inhalation allergens for childhood allergies	DP 3111-1601 E	209
Allercoat ELISA	Mites, animal allergens, grasses, house dust, moulds, tree pollens, herbs and flower pollens	For single determinations of specific IgE against various allergens	EP #####-#####	215
ELISA	Total IgE	ELISA for quantitative determination of human IgE in serum	EV 3840-9601 E	122



To view all EUROIMMUN products for this subject scan the QR code or enter the Quick Link code **q074** at [www.euroimmun.com](http://www.euroimmun.com)

# Specific IgE

Inhalation · Food · Atopy · Insect venoms



For more information on this subject scan the QR code or enter the Quick Link code q003 at [www.euroimmun.com](http://www.euroimmun.com)

## Food

■ **Clinical information:** A food allergy is defined as hypersensitivity to substances which are ingested with the food. This IgE-mediated immune reaction can lead to the corresponding symptoms within a short period after ingesting the food. The foods which most frequently trigger allergic reactions include nuts (particularly peanuts), soya, wheat, fish, milk and eggs. Possible symptoms are burning or itching in the oral cavity, nausea, gastrointestinal spasms, diarrhoea and skin rashes. Severe reactions can also lead to asthma attacks, breathlessness, increased heart rate and to panic attacks and confusion. In rare cases, an anaphylactic shock can occur (e.g. through eating nuts or fish). With a prevalence of 5 to 10%, primary allergic sensitisations to foods play an important role, particularly in babies and infants. Food allergies in adults occur with a prevalence of 1 to 5%.

■ **Diagnostics:** Various **food allergy profiles** are available for the clarification of food allergies. These permit the semi-quantitative and monospecific in vitro determination of human IgE antibodies against the most frequent food allergens in serum or in plasma. Moreover, **country-specific food allergy profiles** are available which have been developed with regard to the regional eating habits.

In addition, the total IgE concentration in the serum can be determined using the Total IgE ELISA.



## Product overview

Method	Substrate	Application	Order number	Page
EUROLINE	Food (f1, f75, f2, f45, f4, f5, f9, f13, f14, f17, f20, f49, f84, f237, f25, f31, f35, f85, f3, f23, CCD)	Efficient screening for IgE antibodies against the most important food allergens	DP 3410-1601 E	210
	Food "Gulf" (f1, f75, f2, f105, f4, f14, f45, fs36, f13, f92, f33, f44, f93, f25, f31, f48, f83, f88, f3, f23, CCD)		DP 3416-1601 E	211
	Food "Turkey 1" (f1, f75, f2, f169, f78, f4, f79, f9, f14, f10, f13, f17, f144, u87, f222, f73, f33, f44, f49, f92, f84, f146, f328, f25, f31, f35, f48, f95, f97, f122, f132, fs14, fs10, fs43, f83, CCD)	Various food profiles with the most important food allergens for particular regions	DP 3420-1601-11 E	211
Allercoat ELISA	Food	For single determinations of specific IgE against various allergens	EP F###-####	220
ELISA	Total IgE	ELISA for quantitative determination of human IgE in serum	EV 3840-9601 E	122



To view all EUROIMMUN products for this subject scan the QR code or enter the Quick Link code **q081** at [www.euroimmun.com](http://www.euroimmun.com)

# Specific IgE

Inhalation · Food · Atopy · Insect venoms



For more information on this subject scan the QR code or enter the Quick Link code q001 at [www.euroimmun.com](http://www.euroimmun.com)

## Atopy

■ **Clinical information:** Atopy is a genetic predisposition to allergic hypersensitivity reactions which can have various clinical manifestations. The allergens responsible for the reaction enter the body either through the air and the mucous membranes (in the case of inhalation allergies), or through food ingestion (in the case of food allergies). However, allergic reactions to foods of plant origin can also be caused by cross-reacting IgE antibodies. These reactions, termed cross-allergies, are based on the structural similarity between proteins which are present in both the food as well as in the corresponding inhalation allergens of plant origin. An example of this is the phenomenon known as "oral allergy syndrome" (OAS). Accordingly, patients with a primary birch pollen allergy can also develop allergic reactions to apple, celery, hazelnut, potato or kiwi. Multiple sensitisations to allergens of different origin are therefore not uncommon.

■ **Diagnostics:** Various **atopy profiles** are available for the clarification of multiple sensitisations. These permit the simultaneous in vitro determination of human IgE antibodies against the most frequent inhalation and food allergens in serum or in plasma. Moreover, **country-specific profiles** are available which take into account the characteristics of the regional allergen exposure.

In addition, the total IgE concentration in the serum can be determined using the Total IgE ELISA.



## Product overview

Method	Substrate	Application	Order number	Page
	Atopy (g6, g12, t3, w6, d1, e1, e2, e3, m2, m6, f1, f2, f3, f4, f9, f14, f17, f31, f35, f49, CCD)	Efficient screening for IgE antibodies against the most important inhalation and food allergens	DP 3710-1601 E	212
EUROLINE	Paediatrics (gx, t3, w6, d1, d2, e1, e2, e3, m2, m3, m6, f1, f75, f2, f3, f76, f77, f78, e204, f4, f9, f14, f13, f17, f31, f35, f49, CCD)	Comprehensive screening profile with the most important inhalation and food allergens for childhood allergies	DP 3712-1601 E	212
Allercoat ELISA	Pharmaceutical drugs, environmental allergens	For single determinations of specific IgE against various allergens	EP #####-#####	215
ELISA	Total IgE	ELISA for quantitative determination of human IgE in serum	EV 3840-9601 E	122



To view all EUROIMMUN products for this subject scan the QR code or enter the Quick Link code q053 at [www.euroimmun.com](http://www.euroimmun.com)

# Specific IgE

Inhalation · Food · Atopy · Insect venoms



For more information on this subject scan the QR code or enter the Quick Link code q028 at [www.euroimmun.com](http://www.euroimmun.com)

## Insect venoms

**Clinical information:** Insect venom allergies are a serious problem, particularly during the summer months. Localised pain, swelling, itching and redness occur as normal reactions to an insect bite. An allergic reaction to an insect bite can occur immediately within minutes or not until hours later. The faster the symptoms occur, the more severe the allergic reaction. With further insect bites, the symptoms can continuously worsen.

According to estimations, around 1 to 5% of the population in Central Europe is sensitised to honey bee and/or wasp venom. In persons who are allergic to insect venoms, severe systemic reactions can occur, which manifest with the formation of urtica, swelling, itching and redness at sites other than the puncture, swelling of the throat and tongue, difficulty in breathing, nausea, gastrointestinal cramps, diarrhoea, neurological deficiencies with confusion, dizziness and gait disorder, as well as raised pulse and fall in blood pressure. An uncontrolled reaction of the immune system to the insect venom can also lead to an anaphylactic shock and thus be life-threatening.



**Diagnostics:** In suspected cases of insect venom allergy, differential diagnostics is recommended for exact identification of the allergy-inducing species. Specific IgE antibodies against bee and wasp venom, as the most common insect venoms, can be detected with the help of a combination test.

Besides the previously used natural **insect venom extracts**, **species-specific allergy components** (defined partial allergen diagnostics, DPA-Dx) are also available for a refined serological diagnosis. In a combined test, both antibodies against the natural bee and wasp venom extracts and against specific molecular antigens for bee and wasp venom (rApi m1, rApi m2, rApi m10, rVes v1 and rVes v5), as well as against CCD (cross-reactive carbohydrate determinant) as a marker for cross-reactivity between the insect venoms can be determined simultaneously in one process. The detection of a primary sensitisation to bee or wasp venom, or of a genuine double sensitisation with comparison to a cross-reaction is made with this test.

In addition, the total IgE concentration in the serum can be determined using the Total IgE ELISA.

## Product overview

Method	Substrate	Application	Order number	Page
EUROLINE	Insect Venoms (i1, i3, CCD)	For reliable diagnosis of insect venom allergies	DP 3720-1601 E	213
	DPA-Dx Insect Venoms 1 (i1, i208, i3, i209, CCD)	Molecular allergy diagnostics with CCD-free major allergens from bee and wasp	DP 3850-1601-1 E	213
	DPA-Dx Insect Venoms 2 (i1, i208, i211, i213, i216, i3, i209, CCD)	Molecular allergy diagnostics with CCD-free major allergens from bee and wasp	DP 3850-1601-2 E	213
Allercoat ELISA	Insects	For single determinations of specific IgE against various allergens	EP I###-0110	229
ELISA	Total IgE	ELISA for quantitative determination of human IgE in serum	EV 3840-9601 E	122



To view all EUROIMMUN products for this subject scan the QR code or enter the Quick Link code **q075** at [www.euroimmun.com](http://www.euroimmun.com)

# Molecular allergy diagnostics

DPA-Dx (Defined partial allergen diagnostics)



For more information on this subject scan the QR code or enter the Quick Link code q039 at [www.euroimmun.com](http://www.euroimmun.com)

## DPA-Dx (Defined partial allergen diagnostics)

■ **Indication:** Most allergy sufferers exhibit specific IgE to several allergen sources. Such multiple sensitisation patterns are often caused by **cross reactions** of IgE antibodies against various inhalation and food allergens.

Pollens are the allergens with the highest prevalence worldwide. Almost 50% of all allergy sufferers have a pollen allergy. As the allergen components of various grasses, trees and herbs are structurally related, pronounced cross-reactions exist with these. Multiple pollen sensitisations are, therefore, not uncommon.


Primary allergic sensitisations play a role in infants and toddlers with a prevalence of 5 to 10%. In infants, mostly allergens from cow's milk, eggs and peanut cause corresponding reactions. A large proportion of the toddlers develop a tolerance against eggs and cow's milk in time. In children, a higher number of allergic reactions to nuts, wheat and fish are observed.

According to estimations, around 1 to 5% of the population in Central Europe is sensitised to honey bee and/or wasp venom. Both venoms possess structurally very similar allergen components, such that a pronounced cross-reactivity also exists here.

■ **Diagnostics:** Conventional in-vitro allergy diagnostic methods for the serological detection of allergen-specific IgE are based on native allergen extracts, which contain an insufficiently defined mixture of different components. The analysis of multiple sensitisations is made difficult by this lack of standardisation of related extracts and underlying cross reactivities are only identified to a limited extent.

The use of **individually purified allergen components** allows the detection of single sensitisations. For the analysis of multiple pollen sensitisations, allergen components of birch (Bet v1, Bet v2, Bet v4 and Bet v6) and grasses (Phlp1, Phlp5, Phlp7 and Phlp12) have been defined as markers for both primary and cross sensitisations. The detection of specific IgE antibodies to these markers therefore enables the **creation of detailed sensitisation profiles** and thus provides important information for the targeted choice of a suitable immunotherapy.

In the area of food allergies, allergen components provide the possibility to estimate the risk of severe reactions and to make statements about the possible development of tolerances. Moreover, primary and pollen-associated food allergies can be differentiated. Therefore, clear statements can be made concerning therapy.



Extract-based diagnostics for insect venom allergies are challenging in that a double sensitisation is diagnosed in vitro in up to 50 % of cases, of which the majority is based on cross-reactions. The largest share of patients (around 80 to 90 %) react in this case to the species-specific major allergens Api m1 and Api m2 of bees or Ves v1 and Ves v5 of wasps. The **combination of allergen extracts and single components** of the respective insect venoms in one test enables the detection of a primary sensitisation to bee or wasp venom or a genuine double sensitisation as opposed to a cross-reaction.

## Product overview

Method	Substrate	Application	Order number	Page
EUROLINE	DPA-Dx Pollen 1 (t3, g6, t215, t216, t220, t225, g205, g215, g210, g212, CCD)	Efficient differential diagnostics of multiple pollen sensitisations	DP 3210-1601-1 E	210
	DPA-Dx Paediatrics (Gal d1, d2, d3, d4, Bos d, d4, d5, d6, d8, Ara h1, h2, h3, h9, Bet v1)	Component-resolved analysis of the most common food allergies in infants	DP 3812-1601-1 E	213
	DPA-Dx Peanut 1 (t215, f422, f423, f424, f427, f429, f444, f445, CCD)	Component-resolved differential diagnostics in peanut allergies	DP 3511-1601-1 E	212
	DPA-Dx Insect Venoms 1 (i1, i208, i3, i209, CCD)	Efficient differential diagnostics with CCD-free major allergens from bee and wasp	DP 3850-1601-1 E	213
	DPA-Dx Insect Venoms 2 (i1, i208, i211, i213, i216, i3, i209, CCD)	Efficient differential diagnostics with CCD-free major allergens from bee and wasp	DP 3850-1601-2 E	213



To view all EUROIMMUN products for this subject scan the QR code or enter the Quick Link code q079 at [www.euroimmun.com](http://www.euroimmun.com)



# **Products for allergy diagnostics**

## Format

- 0501: 5 slides to be incubated with 1 patient serum each
- 0502: 5 slides to be incubated with 2 patient sera each
- 1001: 10 slides to be incubated with 1 patient serum each
- 1601: 16 single test strips
- 0110: 10 paper rings with allergen
- 0125: 25 paper rings with allergen
- 9601: 96 individual break-apart wells (12 microplate strips, 8 wells each)



## Product classification

- DE: Test system EUROASSAY (allergy diagnostics) ..... page 208
- DP: Test system EUROLINE (allergy diagnostics) ..... page 209
- ZP: Allercoat™ 6 System (basic reagents and accessories) ..... page 214
- EP: Allercoat™ 6 System (allergen rings) ..... page 215

For product orders the amount, product code and test name are required. **Test kits** comprise all reagents needed to perform the serological investigation.

Test systems which do not appear in this catalogue can be made to specification. Apart from the customary package sizes and slide formats, special sizes are available as well. Quotations can be provided upon request.

## EUROASSAY for Allergy Diagnostics (Test Systems)

Order No.	Antibodies against	Ig Class	Substrate	Format Slides x Fields
DE 3110-1001 E	inhalation (g1, g3, g6, g12, t2, t3, t4, t7, w1, w8, w6, w9, d1, d2, es4, e2, e3, e6, e1, e82, e84, m1, m2, m3, m6, CCD)	IgE	EUROASSAY strip with allergens	10 x 01
DE 3110-1001-2 E	inhalation 2 (t3, t4, t15, t2, gx2, g12, d1, w6, w1, d2, e1, e2, m6, m2, t7, t11, t23, g6, i6, d70, w21, w9, d71, e3, u85, d72, m3, m1, CCD)	IgE	EUROASSAY strip with allergens	10 x 01
DE 3210-0501-1 E	DPA-Dx pollen 1 (t3, g6, t215, t216, t220, t225, g205, g215, g210, g212, CCD)	IgE	EUROASSAY strip with allergens	05 x 01
DE 3410-1001 E	food (f1, f75, f2, f78, f45, f4, f5, f9, f14, f10, f13, f17, f20, f84, f95, f237, f25, f31, f35, f85, f3, f23, f49, CCD)	IgE	EUROASSAY strip with allergens	10 x 01
DE 3410-1001-2 E	food 2 (f47, f26, f75, f90, f10, f78, f89, u87, f35, f95, f329, f84, f24, f20, f1, f13, f5, f2, f14, f4, f25, f85, f31, f33, f92, f49, f17, f3, CCD)	IgE	EUROASSAY strip with allergens	10 x 01
DE 3712-1001 E	paediatrics/atopy (g6, g12, t3, w6, d1, d2, e1, e2, e3, m2, m3, m6, f1, f75, f3, f2, f76, f77, f78, e204, f4, f9, f14, f13, f17, f31, f35, f49, CCD)	IgE	EUROASSAY strip with allergens	10 x 01
DE 3712-1001-2 E	paediatrics/atopy 2 (t3, t4, t15, t2, gx2, g12, d1, w6, w1, d2, e1, e2, m6, m2, f1, f13, f5, f2, f14, f4, f25, f85, f31, f33, f92, f49, f17, f3, CCD)	IgE	EUROASSAY strip with allergens	10 x 01
DE 3720-0502 E	insect venoms (i1, i3, CCD)	IgE	EUROASSAY strip with allergens	05 x 02
DE 3812-1001-1 E	DPA-Dx paediatrics 1 (t215, f76, f77, f78, f334, e204, f232, f233, f323, f356, f422, f423, f424, f427, CCD)	IgE	EUROASSAY strip with allergens	10 x 01
DE 3850-0501-1 E	DPA-Dx insect venoms 1 (i1, i208, i3, i209, CCD)	IgE	EUROASSAY strip with allergens	05 x 01

**EUROLINE for Allergy Diagnostics (Test Systems)**

Order No.	Antibodies against	Ig Class	Substrate	Format
DP 3110-1601 E	inhalation (g1, g3, g6, g12, t2, t3, t4, t7, w1, w6, w9, d1, d2, e1, e2, e3, m1, m2, m3, m6, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3110-1601-1 E	inhalation 2 (g6, g12, t2, t3, t4, w6, w9, d1, d2, e1, e2, e3, e6, e82, e84, es4, m1, m2, m3, m6, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3111-1601 E	paediatric inhalation (g6, g12, t2, t3, t4, w6, w8, w9, d1, d2, e1, e2, e3, e6, e82, e84, m1, m2, m3, m6, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3112-1601 E	Mediterranean inhalation (g2, g6, t3, t4, t9, t11, t23, t20, w1, w6, w9, w19, d1, d2, d70, e1, e2, e3, m2, m6, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3113-1601 E	inhalation "South East Asia" (ts19, t104, t19, t223, gs1, ds1, i6, u134, e1, e2, es172, e6, e71, e82, e84, ms1, ms4, m5, m12, m45, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3116-1601 E	inhalation "China" (gs23, ts21, t3, t8, t11, t12, t14, t70, ws18, w1, w6, w9, es1, d1, d2, i6, ms5, m1, u73, u80, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3116-1601-2 E	inhalation "China 2" (ds1, h1, i6, e1, e2, ms1, ts20, u80, w1, w6, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3117-1601 E	inhalation "Middle East" (g1, g6, g12, t2, t3, t7, t9, w1, w6, w8, d1, d2, i6, e1, e84, m1, m2, m3, m5, m6, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3118-1601 E	inhalation "Gulf" (g6, g12, t2, t3, t7, t9, w1, w6, d1, d2, i6, e1, e2, e3, e17, m1, m2, m3, m5, m6, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3119-1601-2 E	inhalation "Top Screen" (ds1, e1, e2, e81, gs12, g12, t3, t9, w1, w6, w19, m2, IgE, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3119-1601-3 E	inhalation "Top Screen Plus" (rs1, rs2, d1, d2, e1, e2, e81, g6, g12, t3, t9, w1, w6, w19, m2, IgE, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3119-1601-11 E	inhalation "Turkey 1" (gs12, gs15, gs21, g12, ts23, ts24, t9, t70, ws18, ws19, ws20, d1, d2, i6, es2, es172, e1, e2, e3, e4, e80, e81, e84, ms11, ms12, m1, m2, m3, m6, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3120-1601 E	inhalation "India" (g6, g12, g20, t18, w4, w27, w29, ds1, d2, i6, e1, e2, e11, e85, m3, m37, u81, u126, u129, u140, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3122-1601 E	inhalation "France" (t3, t4, t9, t7, t11, t23, g6, g2, g91, g8, g12, w9, w6, e1, e5, e3, m3, m6, d1, d2, d70, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3122-1601-2 E	inhalation "France 2" (t3, t7, t9, t11, t23, t15, g6, g2, g101, g8, g12, w9, w6, w21, w1, e1, e5, e82, m3, m6, i6, d1, d2, d70, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3123-1601 E	inhalation "Lebanon" (gx2, g2, t9, t23, t11, ts27, ts20, ws19, ws22, w19, w14, d1, d2, e1, e2, es170, i6, u85, m2, m3, m6, ms4, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3123-1601-2 E	inhalation "Lebanon 2" (g6, g2, g12, g101, w7, w19, w21, ws24, t2, t3, t7, t9, t16, t23, ts32, m1, m2, m3, m5, m6, d1, d2, e1, e2, e3, e85, es170, es172, u85, u140, i1, i3, i6, i71, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3124-1601 E	inhalation "Screen South Africa" (ts19, ts29, gs1, ws21, e1, e2, hs12, ms1, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3126-1601 E	inhalation "Mexico" (g2, g5, g6, g10, g14, t3, t7, t19, t210, t102, t14, t15, w1, w4, w6, w8, w10, w11, w14, w15, w100, d1, d2, h1, i6, e1, e2, es4, m1, m2, m3, m6, m20, m5, m11, m14, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01

## EUROLINE for Allergy Diagnostics (Test Systems)

Order No.	Antibodies against	Ig Class	Substrate	Format
DP 3127-1601-1 E	inhalation "Ukraine 1" (d1, d2, e1, e2, e3, e6, e82, e84, t2, t3, t4, t7, g6, g12, w6, w9, m1, m2, m3, m6, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3127-1601-2 E	inhalation "Ukraine 2" (g2, g3, g4, g9, t4, t11, t19, t71, t111, ts31, w7, w8, w11, w12, w14, d72, u85, m1, m2, m3, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3128-1601 E	inhalation trees "Poland" (t1, t2, t3, t4, t5, t7, t15, t12, t14, t16, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3129-1601 E	inhalation grass and weeds "Poland" (g1, g3, g6, g12, w1, w6, w9, w10, w103, w203, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3130-1601 E	inhalation animals "Poland" (e1, e2, e3, e6, e71, e73, e82, e84, es7, es172, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3131-1601 E	inhalation indoor allergens "Poland" (ds1, es2, i6, e7, m1, m2, m3, m5, m6, m37, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3132-1601 E	inhalation "Maghreb" (g2, g3, g6, g8, g12, t3, t5, t7, t9, t11, t15, t18, t19, t23, w1, w4, w6, w7, w9, w10, w21, d1, d2, i6, e1, e2, e3, m3, m5, m6, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3133-1601-1 E	inhalation "Iran 1" (e1, e2, e3, e4, e81, es2, es172, g1, g3, g12, g14, g16, w6, w9, w10, w11, w14, w17, w28, w100, t1, t15, t16, t70, ts22, ts26, i1, i6, h1, d1, d2, m1, m2, m3, m6, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3134-1601-1 E *	tree pollen "South Africa 1" (t2, t3, t7, t8, t9, t11, t12, t14, t15, t16, t18, t19, t23, t70, t111, t112, gs1, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3134-1601-2 E *	tree pollen "South Africa 2" (t2, t3, t7, t8, t9, t11, t12, t14, t15, t16, t18, t19, t23, t70, t111, t112, gs1, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3135-1601-1 E	inhalation "Venezuela 1" (g2, g5, g10, w6, w8, w9, w14, e1, e2, e78, e85, e86, e111, u85, m1, m2, m3, m5, m6, d1, d2, d201, i1, i2, i3, i4, i70, i71, i100, p1, p2, p4, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3210-1601-1 E	DPA-Dx pollen 1 (t3, g6, t215, t216, t220, t225, g205, g215, g210, g212, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3410-1601 E	food (f1, f75, f2, f45, f4, f5, f9, f13, f14, f17, f20, f49, f84, f237, f25, f31, f35, f85, f3, f23, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3410-1601-1 E	food 2 (f1, f75, f2, f78, f4, f5, f14, f10, f13, f17, f20, f49, f84, f95, f25, f31, f35, f85, f3, f23, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3411-1601 E	food "South East Asia 1" (f1, f75, f2, f4, f9, f10, f14, f13, f17, f63, f64, f83, fs10, fs14, f12, f23, f24, f80, f234, f105, f336, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3411-1601-1 E	food "South East Asia 2" (f1, f75, f2, f4, f9, f10, f14, f13, f17, f63, f340, f83, fs10, fs14, f23, f24, f80, f234, f105, f336, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3414-1601 E	food "China" (f1, f2, f4, f7, f27, f88, fs35, f13, f14, fs40, f25, f292, fs42, f23, f234, f3, f41, f56, fs41, fs77, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3414-1601-2 E	food "China 2" (f1, f2, f13, f14, f23, f24, fs33, fs34, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3415-1601 E	food "Middle East" (f1, f75, f2, f78, e204, f4, f14, f45, f13, f17, f20, f33, f49, f92, f25, f31, f85, f48, f88, f89, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01

\*) Currently not available as IVD in the European Union.

## EUROLINE for Allergy Diagnostics (Test Systems)

Order No.	Antibodies against	Ig Class	Substrate	Format
DP 3416-1601 E	food "Gulf" (f1, f75, f2, f105, f4, f14, f45, fs36, f13, f92, f33, f44, f93, f25, f31, f48, f83, f88, f3, f23, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3420-1601-11 E	food "Turkey 1" (f1, f75, f2, f169, f78, f4, f79, f9, f14, f10, f13, f17, f144, u87, f222, f73, f33, f44, f49, f92, f84, f146, f328, f25, f31, f35, f48, f95, f97, f122, f132, fs14, fs10, fs43, f83, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3421-1601 E	food "India" (f2, f75, f168, f4, f9, f14, f13, f36, f49, f50, f35, f38, f48, f244, f83, f89, f74, f240, f23, f24, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3422-1601 E	food "France" (f20, f17, f25, f12, f13, f14, f85, f158, f49, f84, f95, f27, f83, f26, f23, f3, f10, f9, f4, f2, f75, f1, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3422-1601-2 E	food "France 2" (f20, f17, f25, f12, f13, f14, f85, f92, f96, f49, f84, f89, f27, f26, f24, f23, f40, f3, f45, f10, f4, f2, f75, f1, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3423-1601 E	food "Lebanon" (f1, f2, f45, f4, f10, f13, f14, f20, fs39, f12, f49, f95, f33, f44, f86, f31, f35, f25, f63, f83, fs7, fs14, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3423-1601-2 E	food "Lebanon 2" (f1, f2, f23, f24, f40, f41, fs79, f4, f9, f10, f12, f13, f14, f20, fs39, f45, f27, f83, f44, f49, f84, f95, fs80, f25, f31, f35, f47, f48, f86, f96, f159, f292, f273, fs78, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3424-1601 E	food "Screen South Africa" (f1, f2, f4, f13, f14, fs36, f3, fs12, fs58, fs53, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3425-1601 E	food "Mexico" (f1, f75, f2, f78, f95, f96, f44, fs32, f4, f7, f9, f45, f13, f14, f20, fs35, f12, f15, f49, f292, f25, f31, f35, f92, f216, f191, f263, f105, f284, f83, f26, f27, f24, f40, fs12, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3426-1601-1 E	food "Ukraine 1" (f1, f75, f2, f78, f4, f5, f14, f10, f13, f256, f17, f20, f49, f95, f25, f31, f35, f85, f3, f24, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3426-1601-2 E	food "Ukraine 2" (f4, f6, f9, f11, f13, f361, f45, f47, f73, f81, fs32, f25, f35, f26, f27, f83, f3, f40, f41, f206, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3427-1601 E	food dairy products and nuts "Poland" (f1, f2, f75, f78, f13, f256, f17, f20, f73, f336, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3428-1601 E	food flour and meat "Poland" (f4, f5, f7, f9, f26, f27, f83, f79, f3, f24, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3429-1601 E	food vegetables "Poland" (f10, f14, f86, f25, f31, f35, f85, f46, f244, f292, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3430-1601 E	food fruits "Poland" (f44, f49, f84, f92, f95, f97, f122, f237, f329, fs32, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3431-1601 E	food "Maghreb" (f1, f75, f2, f78, f4, f9, f10, f12, f13, f14, f17, f20, f33, f44, f49, f92, f25, f35, f47, f48, f85, f45, f73, f27, f83, f3, f40, f41, f23, f24, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3432-1601-1 E	food "Iran 1" (f1, f75, f2, f78, f4, f5, f9, fs13, f13, f14, f17, f20, f144, f256, f79, f44, f49, f50, f84, f92, f93, f95, f97, f87, fs32, f25, f35, f46, f47, f85, f262, f83, fs28, f24, fs12, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3433-1601-1 E	food "Venezuela 1" (f1, f2, f3, f4, f5, f6, f7, f9, f13, f14, f256, f17, f20, f24, f25, f26, f27, f92, f31, f32, f33, f35, f37, f40, f41, f43, f44, f45, f49, f233, f72, f75, f76, f77, f78, f79, f81, f83, f105, f160, f205, f219, f292, f232, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3510-1601-1 E	DPA-Dx milk 1 (f2, f76, f77, f78, f334, e204, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01

## EUROLINE for Allergy Diagnostics (Test Systems)

Order No.	Antibodies against	Ig Class	Substrate	Format
DP 3511-1601-1 E	DPA-Dx peanuts 1 (t215, f422, f423, f424, f429, f445, f444, f427, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3512-1601-1 E	DPA-Dx food 1 (f1, f2, f13, t215, f76, f77, f78, f334, e204, f232, f233, f323, f356, f422, f423, f424, f427, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3701-1601 E	atopy "Top Screen" (rs1, rs2, fx5, fs52, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3702-1601-1 E	atopy "Chile 1" (d1, d2, e1, e2, es2, g2, g4, g5, m2, m3, m6, t11, t14, t15, t19, w9, w100, f1, f75, f78, f2, f40, f41, fs12, f4, f9, f13, f14, f25, f33, f49, f92, f83, f96, f105, f292, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3702-1601-2 E	atopy "Chile 2" (d1, d2, e1, e2, es2, gs2, m5, i6, t1, t7, t11, t14, t19, ts30, ws23, f1, f2, f4, f6, f7, f13, f17, f36, f256, i14, f74, f105, f292, f26, f27, f83, f88, f3, f40, f41, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3703-1601-1 E	atopy "Ukraine 1" (d1, d2, e1, e2, IgE, h1, is6, m6, g6, g12, ts3, t7, w1, w6, f1, f2, f14, f95, f23, fs59, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3703-1601-2 E	atopy "Ukraine 2" (d1, d2, e1, e2, e3, g6, t3, w6, m2, m3, f1, f75, f2, f3, f4, f13, f14, f14, f31, f35, f49, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3704-1601-1 E	atopy "Venezuela 1" (f1, f75, f3, f24, f40, f308, f4, f5, f6, f7, f9, f14, f292, f79, f13, f17, f20, f2, f76, f77, f78, f105, f218, f25, f32, f33, f44, f49, f72, f26, f27, f83, f155, f45, u85, d1, d2, d201, e1, e2, e85, m1, m2, m3, m5, m6, g2, w8, i1, i3, i70, i71, i100, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3705-1601-1 E	atopy "Korea 1" (g1, g2, g3, g6, g7, g12, e1, e2, t4, t7, t12, t16, t19, t71, ts3, ts11, ts14, ts15, u80, f1, f2, f23, f24, f206, f14, f95, w1, w6, w7, w8, w11, w12, w14, i6, d1, d2, h1, m1, m2, m3, m6, IgE, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3705-1601-2 E	atopy "Korea 2" (f1, f2, f81, f26, f27, f83, f3, f40, f41, f206, f23, f24, f4, f6, f9, f11, f13, f14, fs15, f95, f25, f47, f48, f45, d1, d2, d70, e1, e2, h1, m2, m3, m5, m6, i6, u80, g12, ts3, t7, w1, w6, IgE, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3706-1601-1 E	atopy "Caribbean 1" (g6, g12, t3, t20, w6, e1, e2, e3, d1, m2, m6, f1, f2, f79, f4, f9, f14, f13, f31, f35, f292, f49, f362, f64, f3, f24, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3707-1601-1 E	atopy "Indonesia 1" (g2, g6, gs1, t19, t104, t223, d1, d2, d4, d72, d73, d201, e1, e2, e3, e204, es2, u134, i1, i6, m5, ms1, f4, f79, f13, f14, f17, f20, f45, f336, f1, f2, f76, f77, f78, fs10, f3, f40, f41, f23, f24, f80, f58, f63, f64, f83, f88, f81, f25, f47, f44, f84, f74, f105, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3710-1601 E	atopy (g6, g12, t3, w6, d1, e1, e2, e3, m2, m6, f1, f2, f3, f4, f9, f14, f17, f31, f35, f49, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3710-1601-2 E	atopy 3 (g6, t3, t4, w6, d1, d2, e1, e2, e3, m2, m3, f1, f75, f2, f3, f4, f13, f14, f31, f49, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3711-1601 E	pollen-food cross reactions (g6, t3, w6, f4, f5, f13, f17, f20, f48, f89, f271, f275, f44, f49, f348, f237, f328, f31, f35, f85, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3712-1601 E	paediatrics (gx, t3, w6, d1, d2, e1, e2, e3, m2, m3, m6, f1, f75, f2, f3, f76, f77, f78, e204, f4, f9, f14, f13, f17, f31, f35, f49, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3713-1601 E	atopy "China"	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3713-6401 E	(ts20, w1, w6, ds1, h1, e1, e2, i6, ms1, u80, f1, f2, f13, f14, f27, f88, fs33, fs34, f24, f23, CCD)			64 x 01
DP 3713-1601-3 E	atopy "China 3" (ts22, w6, us1, ds1, es1, h1, ms5, f245, fs9, fs43, fs56, fs34, fs44, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01

**EUROLINE for Allergy Diagnostics (Test Systems)**

Order No.	Antibodies against	Ig Class	Substrate	Format
DP 3713-1601-4 E	atopy "China 4" (ds1, h1, i6, e1, e2, ms1, ts20, u80, w1, w6, f1, f2, f13, fs33, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3713-1601-5 E	atopy "China 5" (ds1, h1, i6, e1, e2, ms1, w1, w6, f1, f2, f13, f14, f23, f24, fs33, u80, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3713-1601-7 E	atopy "China 7" (ds1, h1, i6, e1, e2, ms1, w1, w6, f1, f2, fs33, u80, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3713-1601-8 E DP 3713-6401-8 E	atopy "China 8" (IgE, ts20, w1, w6, ds1, h1, e1, e2, i6, ms1, u80, f1, f2, f13, f14, f27, f88, fs33, f24, f23, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01 64 x 01
DP 3716-1601 E	atopy "Peru" (e1, e2, e3, e7, e11, e82, e84, es172, i1, i3, i6, m1, m2, m3, m6, t18, u85, w29, w28, d1, d2, f75, f79, f1, f2, f23, f206, f13, f14, f104, f105, f26, fs32, f44, f49, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3717-1601 E	"Mix France" (f85, f14, f13, f12, f4, f23, f3, f2, f1, f17, f84, f95, t3, t9, t7, t11, t23, g6, w9, w1, w6, w19, e1, e5, m3, m6, d1, d2, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3717-1601-2 E	"Mix France 2" (f85, f14, f13, f12, f4, f24, f3, f2, f1, f17, f49, f84, t3, t7, t9, t15, t23, g6, w9, w6, w21, w1, e1, e5, m3, m6, d1, d2, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3718-1601 E	atopy "Lithuania" (w6, w9, w103, w203, t3, t4, t11, g6, gs21, e1, e2, e3, e4, e82, es2, ds1, ms1, m5, f73, f245, f2, f81, f13, f14, f17, f256, f4, f12, f15, f44, f32, f26, f27, f83, f49, f3, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3720-1601 E	insect venoms (i1, i3, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3790-1601 E	atopy screen (d1, d2, i1, i3, i6, h1, e1, e2, e3, m1, m2, m3, m6, g1, g3 g6, g12, t2, t3, t4, t7, t23, w1, w6, w9, u85, f25, f31, f35, f85, f1, f75, f2, f3, f23, f24, e204, f76, f77, f78, f27, f88, f45, f4, f5, f9, f14, f10, f13, f17, f20, f49, f84, f237, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3812-1601-1 E	DPA-Dx paediatrics 1 (t215, f76, f77, f78, f334, e204, f232, f233, f323, f356, f422, f423, f424, f427, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3850-1601-1 E	DPA-Dx insect venoms 1 (i1, i208, i3, i209, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01
DP 3850-1601-2 E	DPA-Dx insect venoms 2 (i1, i208, i213, i216, i3, i209, i211, CCD)	IgE	membrane strip with allergens (EUROLINE)	16 x 01

## Allercoat™ 6 Systems: Basic Reagents and Accessories

Order No.	Previous Order No.	Reagent	Format
ZP 3000-0100	22000	conjugate pack	for 100 analyses
ZP 3001-0060	22001	reference set	for 5 reference curves in double determination
ZP 3002-0511	22002	washing solution (concentrate)	for 5 x 500 ml
ZP 3003-0511	22003	substrate buffer (incl. 5 substrate tablets)	5 x 13.5 ml, ready for use
ZP 3004-0511	22004	stop solution	5 x 11 ml, ready for use
ZP 3008-0125	22008	positive control serum (G5)	2.5 ml, ready for use
ZP 3009-0125	22009	negative control serum (G5)	2.5 ml, ready for use
ZP 3020-48001		microplate (only in combination with further reagents)	5 microplates with 96 wells each

**Allercoat™ 6 System: Allergen Rings (Pharmaceutical Drugs)**

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP C000-0110	22414	HSA: Human serum albumin	paper rings with allergen	01 x 10
EP C001-0110	22250	C1: Penicilloyl G - HSA	paper rings with allergen	01 x 10
EP C002-0110	22251	C2: Penicilloyl V - HSA	paper rings with allergen	01 x 10
EP C003-0110	22252	C3: A.C.T.H.	paper rings with allergen	01 x 10
EP C004-0110	22253	C4: Chymopapain (Disease)	paper rings with allergen	01 x 10
EP C005-0110	22254	C5: Chymopapain	paper rings with allergen	01 x 10
EP C006-0110	22255	C6: Cephalosporin - HSA	paper rings with allergen	01 x 10
EP C007-0110	22256	C7: Bactrim - HSA	paper rings with allergen	01 x 10
EP C025-0110	22570	C25: Allopurinol - HSA	paper rings with allergen	01 x 10
EP C049-0110	22575	C49: Propicillin - HSA	paper rings with allergen	01 x 10
EP C054-0110	22576	C54: Cefalotin - HSA	paper rings with allergen	01 x 10
EP C060-0110	22577	C60: Gentamycin - HSA	paper rings with allergen	01 x 10
EP C061-0110	22578	C61: Erythromycin - HSA	paper rings with allergen	01 x 10
EP C062-0110	22579	C62: Doxycyclin - HSA	paper rings with allergen	01 x 10
EP C064-0110	22580	C64: Piperacillin - HSA	paper rings with allergen	01 x 10
EP C065-0110	22581	C65: Phenylbutazone - HSA	paper rings with allergen	01 x 10
EP C067-0110	22582	C67: Cloxacillin - HSA	paper rings with allergen	01 x 10
EP C068-0110	22583	C68: Articain - HSA	paper rings with allergen	01 x 10
EP C073-0110	22270	C73: Human insulin	paper rings with allergen	01 x 10
EP C075-0110	22584	C75: Bromelain - HSA	paper rings with allergen	01 x 10
EP C077-0110	22585	C77: Piroxicam - HSA	paper rings with allergen	01 x 10
EP C078-0110	22586	C78: Ibuprofen - HSA	paper rings with allergen	01 x 10
EP C079-0110	22587	C79: Diclofenac - HSA	paper rings with allergen	01 x 10
EP C080-0110	22588	C80: Tetanus	paper rings with allergen	01 x 10
EP C081-0110	22589	C81: Theophylline/Aminophylline - HSA	paper rings with allergen	01 x 10
EP C082-0110	22590	C82: Lidocaine - HSA	paper rings with allergen	01 x 10
EP C083-0110	22553	C83: Novocaine - HSA	paper rings with allergen	01 x 10
EP C086-0110	22591	C86: Benzocaine - HSA	paper rings with allergen	01 x 10
EP C088-0110	22554	C88: Mepivacaine - HSA	paper rings with allergen	01 x 10
EP C089-0110	22592	C89: Bupivacaine - HSA	paper rings with allergen	01 x 10
EP C090-0110	22593	C90: Propyphenazone - HSA	paper rings with allergen	01 x 10
EP C091-0110	22594	C91: Metamizol - HSA	paper rings with allergen	01 x 10
EP C093-0110	22595	C93: Indomethacin - HSA	paper rings with allergen	01 x 10
EP C094-0110	22555	C94: Tobramycin - HSA	paper rings with allergen	01 x 10
EP C095-0110	22556	C95: Neomycin - HSA	paper rings with allergen	01 x 10
EP C096-0110	22596	C96: Ambroxol - HSA	paper rings with allergen	01 x 10
EP C097-0110	22597	C97: Bromhexine - HSA	paper rings with allergen	01 x 10
EP C098-0110	22598	C98: Codeine - HSA	paper rings with allergen	01 x 10
EP C099-0110	22599	C99: L-Tyroxin - HSA	paper rings with allergen	01 x 10
EP C100-0110	22559	C100: Prilocain - HSA	paper rings with allergen	01 x 10

## Allercoat™ 6 System: Allergen Rings (Pharmaceutical Drugs)

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP C101-0110	22560	C101: Phenobarbital - HSA	paper rings with allergen	01 x 10
EP C102-0110	22549	C102: Thiopental - HSA	paper rings with allergen	01 x 10
EP C103-0110	22550	C103: Isoprenaline - HSA	paper rings with allergen	01 x 10
EP C104-0110	22561	C104: Clindamycin - HSA	paper rings with allergen	01 x 10
EP C106-0110	22562	C106: Vitamin B1 - HSA	paper rings with allergen	01 x 10
EP C107-0110	22563	C107: Captopril - HSA	paper rings with allergen	01 x 10
EP C108-0110	22551	C108: Ciprofloxacin - HSA	paper rings with allergen	01 x 10
EP C109-0110	22564	C109: Vitamin B6 - HSA	paper rings with allergen	01 x 10
EP C110-0110	22566	C110: Naproxen - HSA	paper rings with allergen	01 x 10
EP C111-0110	22567	C111: Phenacetin - HSA	paper rings with allergen	01 x 10
EP C112-0110	22568	C112: Diphtheria toxin	paper rings with allergen	01 x 10
EP C200-0110	22257	C200: Paracetamol - HSA	paper rings with allergen	01 x 10
EP C201-0110	22258	C201: Acetylsalicylic acid - HSA	paper rings with allergen	01 x 10
EP C203-0110	22259	C203: Ampicillin - HSA	paper rings with allergen	01 x 10
EP C204-0110	22260	C204: Amoxicillin - HSA	paper rings with allergen	01 x 10
EP C205-0110	22261	C205: Tetracycline - HSA	paper rings with allergen	01 x 10
EP C206-0110	22262	C206: Streptomycin - HSA	paper rings with allergen	01 x 10
EP C207-0110	22263	C207: Trimethoprim - HSA	paper rings with allergen	01 x 10
EP C208-0110	22264	C208: Sulphonmethoxazol - HSA	paper rings with allergen	01 x 10
EP C209-0110	22265	C209: Pyrazolone - HSA	paper rings with allergen	01 x 10
EP C210-0110	22266	C210: Muscle relaxant - HSA	paper rings with allergen	01 x 10
EP C223-0110	22887	C223: Benzodiazepine - HSA	paper rings with allergen	01 x 10
EP C300-0110	22267	C300: Procaine - HSA	paper rings with allergen	01 x 10
EP C320-0110	22572	C320: Nifurantan - HSA	paper rings with allergen	01 x 10
EP C379-0110	22888	C379: Butyrophenone/Haloperidol - HSA	paper rings with allergen	01 x 10
EP C424-0110	22573	C424: Prednisolone - HSA	paper rings with allergen	01 x 10
EP C430-0110	22574	C430: Pyrimethamine - HSA	paper rings with allergen	01 x 10
EP C431-0110	22565	C431: Ascorbic acid - HSA	paper rings with allergen	01 x 10

**Allercoat™ 6 System: Allergen Rings (Mites)**

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP D001-0125	22050	D1: <i>Dermatophagoides pteronyssinus</i>	paper rings with allergen	01 x 25
EP D002-0125	22051	D2: <i>Dermatophagoides farinae</i>	paper rings with allergen	01 x 25
EP D004-0110	22600	D4: <i>Dermatophagoides microceras</i>	paper rings with allergen	01 x 10
EP D070-0110	22271	D70: <i>Acarus siro</i> (storage mite)	paper rings with allergen	01 x 10
EP D071-0110	22272	D71: <i>Glycyphagus destructor</i>	paper rings with allergen	01 x 10
EP D072-0110	22273	D72: <i>Tyrophagus putrescentiae</i>	paper rings with allergen	01 x 10
EP D073-0110	22274	D73: <i>Glycyphagus domesticus</i>	paper rings with allergen	01 x 10
EP D074-0110	22275	D74: <i>Euroglyphus maynei</i>	paper rings with allergen	01 x 10

<b>Allercoat™ 6 System: Allergen Rings (Animal Allergens)</b>				
<b>Order No.</b>	<b>Previous Order No.</b>	<b>Antibodies of Class IgE against</b>	<b>Substrate</b>	<b>Format</b>
EP E001-0125	22052	<b>E1: Cat epithelium</b>	paper rings with allergen	01 x 25
EP E002-0125	22053	<b>E2: Dog hair + epithelium</b>	paper rings with allergen	01 x 25
EP E003-0125	22054	<b>E3: Horse epithelium</b>	paper rings with allergen	01 x 25
EP E004-0110	22276	<b>E4: Cow epithelium</b>	paper rings with allergen	01 x 10
EP E005-0125	22055	<b>E5: Dog squames</b>	paper rings with allergen	01 x 25
EP E006-0125	22056	<b>E6: Guinea pig epithelium</b>	paper rings with allergen	01 x 25
EP E007-0110	22277	<b>E7: Pigeon droppings</b>	paper rings with allergen	01 x 10
EP E008-0110	22557	<b>E8: Mink fur and epithelium</b>	paper rings with allergen	01 x 10
EP E011-0110	22722	<b>E11: Pigeon feathers</b>	paper rings with allergen	01 x 10
EP E012-0110	22725	<b>E12: Pigeon protein</b>	paper rings with allergen	01 x 10
EP E013-0110	22726	<b>E13: Pigeon serum</b>	paper rings with allergen	01 x 10
EP E014-0110	22603	<b>E14: Canary serum</b>	paper rings with allergen	01 x 10
EP E015-0110	22727	<b>E15: Chicken serum</b>	paper rings with allergen	01 x 10
EP E016-0110	22604	<b>E16: Parrot serum</b>	paper rings with allergen	01 x 10
EP E017-0110	22728	<b>E17: Camel hair (native) + epithelium</b>	paper rings with allergen	01 x 10
EP E018-0110	22605	<b>E18: Canary droppings</b>	paper rings with allergen	01 x 10
EP E019-0110	22729	<b>E19: Goose droppings</b>	paper rings with allergen	01 x 10
EP E020-0110	22606	<b>E20: Chicken droppings</b>	paper rings with allergen	01 x 10
EP E030-0110	22730	<b>E30: Human hair</b>	paper rings with allergen	01 x 10
EP E050-0110	22731	<b>E50: Finch feathers + droppings</b>	paper rings with allergen	01 x 10
EP E052-0110	22607	<b>E52: Rabbit/hare epithelium</b>	paper rings with allergen	01 x 10
EP E070-0125	22057	<b>E70: Goose feathers</b>	paper rings with allergen	01 x 25
EP E071-0110	22280	<b>E71: Mouse epithelium</b>	paper rings with allergen	01 x 10
EP E072-0110	22281	<b>E72: Mouse urine proteins</b>	paper rings with allergen	01 x 10
EP E073-0110	22282	<b>E73: Rat epithelium</b>	paper rings with allergen	01 x 10
EP E074-0110	22283	<b>E74: Rat urine proteins</b>	paper rings with allergen	01 x 10
EP E075-0110	22284	<b>E75: Rat serum proteins</b>	paper rings with allergen	01 x 10
EP E076-0110	22285	<b>E76: Mouse serum proteins</b>	paper rings with allergen	01 x 10
EP E077-0110	22286	<b>E77: Budgerigar droppings</b>	paper rings with allergen	01 x 10
EP E078-0125	22058	<b>E78: Budgerigar feathers</b>	paper rings with allergen	01 x 25
EP E079-0110	22287	<b>E79: Budgerigar serum proteins</b>	paper rings with allergen	01 x 10
EP E080-0110	22288	<b>E80: Goat epithelium</b>	paper rings with allergen	01 x 10
EP E081-0110	22289	<b>E81: Sheep epithelium</b>	paper rings with allergen	01 x 10
EP E082-0110	22290	<b>E82: Rabbit epithelium</b>	paper rings with allergen	01 x 10
EP E083-0110	22291	<b>E83: Pig epithelium</b>	paper rings with allergen	01 x 10
EP E084-0110	22292	<b>E84: Hamster epithelium</b>	paper rings with allergen	01 x 10
EP E085-0110	22293	<b>E85: Chicken feathers</b>	paper rings with allergen	01 x 10
EP E086-0125	22059	<b>E86: Duck feathers</b>	paper rings with allergen	01 x 25
EP E089-0110	22732	<b>E89: Mouse droppings</b>	paper rings with allergen	01 x 10
EP E090-0110	22294	<b>E90: Canary feathers</b>	paper rings with allergen	01 x 10

**Allercoat™ 6 System: Allergen Rings (Animal Allergens)**

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP E091-0110	22295	<b>E91: Parrot feathers</b>	paper rings with allergen	01 x 10
EP E093-0110	22608	<b>E93: Love bird serum</b>	paper rings with allergen	01 x 10
EP E094-0110	22737	<b>E94: Ring-necked parakeet droppings</b>	paper rings with allergen	01 x 10
EP E096-0110	22610	<b>E96: Love bird droppings</b>	paper rings with allergen	01 x 10
EP E097-0110	22609	<b>E97: Parrot droppings</b>	paper rings with allergen	01 x 10
EP E098-0110	22738	<b>E98: Chinchilla epithelium + hair</b>	paper rings with allergen	01 x 10
EP E099-0110	22739	<b>E99: Rabbit droppings</b>	paper rings with allergen	01 x 10
EP E100-0110	22558	<b>E100: Duck droppings</b>	paper rings with allergen	01 x 10
EP E103-0110	22721	<b>E103: Wild-pig epithelium + hair</b>	paper rings with allergen	01 x 10
EP E104-0110	22601	<b>E104: Horse hair</b>	paper rings with allergen	01 x 10
EP E110-0110	22723	<b>E110: Ring-necked parakeet feathers</b>	paper rings with allergen	01 x 10
EP E111-0110	22602	<b>E111: Goose down</b>	paper rings with allergen	01 x 10
EP E113-0110	22724	<b>E113: Zebra finch feathers</b>	paper rings with allergen	01 x 10
EP E120-0110	22733	<b>E120: Rat droppings</b>	paper rings with allergen	01 x 10
EP E121-0110	22736	<b>E121: Turkey feathers</b>	paper rings with allergen	01 x 10
EP E204-0110	22278	<b>E204: Bovine serum albumin</b>	paper rings with allergen	01 x 10
EP E205-0110	22279	<b>E205: Horse serum albumin</b>	paper rings with allergen	01 x 10
EP ES11-0125	22150	<b>ES11: Animals (E1, E3, E4, E5)</b>	paper rings with allergen mixture	01 x 25
EP ES170-0110	22950	<b>ES170: Rodents (E6, E71, E73, E82, E84)</b>	paper rings with allergen mixture	01 x 10
EP ES171-0125	22151	<b>ES171: Feathers (E70, E85, E86)</b>	paper rings with allergen mixture	01 x 25
EP ES172-0110	22963	<b>ES172: Cage birds (E78, E90, E91, E110, E 113)</b>	paper rings with allergen mixture	01 x 10
EP ES173-0125	22171	<b>ES173: Bed feathers (E86, E70)</b>	paper rings with allergen mixture	01 x 25

<b>Allercoat 6™-System: Allergen Rings (Food)</b>				
<b>Order No.</b>	<b>Previous Order No.</b>	<b>Antibodies of Class IgE against</b>	<b>Substrate</b>	<b>Format</b>
EP F001-0125	22060	<b>F1: Hens egg white</b>	paper rings with allergen	01 x 25
EP F002-0125	22061	<b>F2: Cows milk protein</b>	paper rings with allergen	01 x 25
EP F003-0125	22062	<b>F3: Codfish</b>	paper rings with allergen	01 x 25
EP F004-0125	22063	<b>F4: Wheat flour</b>	paper rings with allergen	01 x 25
EP F005-0125	22064	<b>F5: Rye flour</b>	paper rings with allergen	01 x 25
EP F006-0110	22296	<b>F6: Barley flour</b>	paper rings with allergen	01 x 10
EP F007-0110	22297	<b>F7: Oat flour</b>	paper rings with allergen	01 x 10
EP F008-0110	22300	<b>F8: Corn flour (Maize flour)</b>	paper rings with allergen	01 x 10
EP F009-0110	22301	<b>F9: Rice</b>	paper rings with allergen	01 x 10
EP F010-0110	22302	<b>F10: Sesame seeds, ground</b>	paper rings with allergen	01 x 10
EP F011-0110	22311	<b>F11: Buckwheat flour</b>	paper rings with allergen	01 x 10
EP F012-0110	22314	<b>F12: Pea</b>	paper rings with allergen	01 x 10
EP F013-0125	22066	<b>F13: Peanuts</b>	paper rings with allergen	01 x 25
EP F014-0125	22067	<b>F14: Soybean</b>	paper rings with allergen	01 x 25
EP F015-0110	22315	<b>F15: White beans</b>	paper rings with allergen	01 x 10
EP F016-0110	22316	<b>F16: Walnut</b>	paper rings with allergen	01 x 10
EP F017-0125	22068	<b>F17: Hazelnut</b>	paper rings with allergen	01 x 25
EP F018-0110	22317	<b>F18: Brazil nut</b>	paper rings with allergen	01 x 10
EP F019-0110	22318	<b>F19: Sweet chestnut</b>	paper rings with allergen	01 x 10
EP F020-0110	22319	<b>F20: Almond</b>	paper rings with allergen	01 x 10
EP F021-0110	22321	<b>F21: Herring</b>	paper rings with allergen	01 x 10
EP F022-0110	22323	<b>F22: Trout</b>	paper rings with allergen	01 x 10
EP F023-0110	22324	<b>F23: Crab</b>	paper rings with allergen	01 x 10
EP F024-0110	22330	<b>F24: Prawn, North Sea</b>	paper rings with allergen	01 x 10
EP F025-0125	22069	<b>F25: Tomato</b>	paper rings with allergen	01 x 25
EP F026-0125	22070	<b>F26: Pork, raw</b>	paper rings with allergen	01 x 25
EP F027-0125	22071	<b>F27: Beef, raw</b>	paper rings with allergen	01 x 25
EP F029-0110	22350	<b>F29: Banana</b>	paper rings with allergen	01 x 10
EP F030-0110	22353	<b>F30: Grapefruit</b>	paper rings with allergen	01 x 10
EP F031-0125	22072	<b>F31: Carrot</b>	paper rings with allergen	01 x 25
EP F032-0110	22354	<b>F32: Lemon</b>	paper rings with allergen	01 x 10
EP F033-0125	22073	<b>F33: Orange</b>	paper rings with allergen	01 x 25
EP F034-0110	22646	<b>F34: Tangerine</b>	paper rings with allergen	01 x 10
EP F035-0125	22074	<b>F35: Potato</b>	paper rings with allergen	01 x 25
EP F036-0110	22355	<b>F36: Coconut</b>	paper rings with allergen	01 x 10
EP F037-0110	22356	<b>F37: Mussel</b>	paper rings with allergen	01 x 10
EP F038-0110	22647	<b>F38: Spinach</b>	paper rings with allergen	01 x 10
EP F039-0110	22648	<b>F39: Cabbage, white</b>	paper rings with allergen	01 x 10
EP F040-0110	22357	<b>F40: Tuna</b>	paper rings with allergen	01 x 10
EP F041-0110	22358	<b>F41: Salmon</b>	paper rings with allergen	01 x 10

**Allercoat 6™-System: Allergen Rings (Food)**

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP F043-0110	22359	F43: Brewer's yeast <i>Saccharomyces carlsbergensis</i>	paper rings with allergen	01 x 10
EP F044-0110	22360	F44: Strawberry	paper rings with allergen	01 x 10
EP F045-0110	22361	F45: Baker's yeast <i>Saccharomyces cerevisiae</i>	paper rings with allergen	01 x 10
EP F046-0110	22649	F46: Pepper, green	paper rings with allergen	01 x 10
EP F047-0110	22362	F47: Garlic	paper rings with allergen	01 x 10
EP F049-0110	22363	F49: Apple	paper rings with allergen	01 x 10
EP F050-0110	22650	F50: Grape	paper rings with allergen	01 x 10
EP F051-0110	22651	F51: Soy	paper rings with allergen	01 x 10
EP F055-0110	22364	F55: Eel	paper rings with allergen	01 x 10
EP F056-0110	22365	F56: Sea perch	paper rings with allergen	01 x 10
EP F057-0110	22366	F57: Goose meat	paper rings with allergen	01 x 10
EP F058-0110	22367	F58: Duck meat	paper rings with allergen	01 x 10
EP F062-0110	22652	F62: Cauliflower, cooked	paper rings with allergen	01 x 10
EP F063-0110	22653	F63: Beef, cooked	paper rings with allergen	01 x 10
EP F066-0110	22812	F66: Leek	paper rings with allergen	01 x 10
EP F068-0110	22654	F68: Ovomucoid	paper rings with allergen	01 x 10
EP F070-0110	22655	F70: Swiss cheese	paper rings with allergen	01 x 10
EP F071-0110	22368	F71: Spiny lobster	paper rings with allergen	01 x 10
EP F072-0110	22369	F72: Pineapple	paper rings with allergen	01 x 10
EP F073-0110	22370	F73: Cocoa	paper rings with allergen	01 x 10
EP F074-0110	22371	F74: Coffee	paper rings with allergen	01 x 10
EP F075-0125	22075	F75: Hens egg yolk	paper rings with allergen	01 x 25
EP F076-0110	22372	F76: alpha-lactalbumin	paper rings with allergen	01 x 10
EP F077-0110	22373	F77: beta-lactoglobulin	paper rings with allergen	01 x 10
EP F078-0110	22374	F78: Casein	paper rings with allergen	01 x 10
EP F079-0110	22375	F79: Gluten (Gliadin)	paper rings with allergen	01 x 10
EP F080-0110	22376	F80: Lobster	paper rings with allergen	01 x 10
EP F081-0110	22377	F81: Cheddar cheese	paper rings with allergen	01 x 10
EP F082-0110	22378	F82: Blue cheese	paper rings with allergen	01 x 10
EP F083-0110	22379	F83: Chicken	paper rings with allergen	01 x 10
EP F084-0110	22380	F84: Kiwi fruit	paper rings with allergen	01 x 10
EP F085-0110	22381	F85: Celery	paper rings with allergen	01 x 10
EP F086-0110	22382	F86: Parsley	paper rings with allergen	01 x 10
EP F087-0110	22383	F87: Melon	paper rings with allergen	01 x 10
EP F088-0110	22384	F88: Mutton / Lamb	paper rings with allergen	01 x 10
EP F089-0110	22385	F89: Mustard	paper rings with allergen	01 x 10
EP F090-0110	22386	F90: Pepper, black	paper rings with allergen	01 x 10
EP F091-0110	22387	F91: Paprika	paper rings with allergen	01 x 10
EP F092-0110	22388	F92: Nutmeg	paper rings with allergen	01 x 10
EP F093-0110	22389	F93: Mango	paper rings with allergen	01 x 10

<b>Allercoat 6™-System: Allergen Rings (Food)</b>				
<b>Order No.</b>	<b>Previous Order No.</b>	<b>Antibodies of Class IgE against</b>	<b>Substrate</b>	<b>Format</b>
EP F094-0110	22390	<b>F94: Pear</b>	paper rings with allergen	01 x 10
EP F095-0110	22391	<b>F95: Peach</b>	paper rings with allergen	01 x 10
EP F096-0110	22392	<b>F96: Avocado</b>	paper rings with allergen	01 x 10
EP F097-0110	22393	<b>F97: Cherry</b>	paper rings with allergen	01 x 10
EP F098-0110	22394	<b>F98: Millet</b>	paper rings with allergen	01 x 10
EP F099-0110	22395	<b>F99: Spelt</b>	paper rings with allergen	01 x 10
EP F100-0110	22303	<b>F100: Lettuce</b>	paper rings with allergen	01 x 10
EP F101-0110	22304	<b>F101: Amaranth seeds</b>	paper rings with allergen	01 x 10
EP F102-0110	22305	<b>F102: Malt, made of barley</b>	paper rings with allergen	01 x 10
EP F103-0110	22306	<b>F103: Pecan nut</b>	paper rings with allergen	01 x 10
EP F104-0110	22307	<b>F104: Quinoa</b>	paper rings with allergen	01 x 10
EP F105-0125	22065	<b>F105: Chocolate</b>	paper rings with allergen	01 x 25
EP F106-0110	22308	<b>F106: Kohlrabi</b>	paper rings with allergen	01 x 10
EP F108-0110	22309	<b>F108: Onion</b>	paper rings with allergen	01 x 10
EP F109-0110	22310	<b>F109: Cotton seeds</b>	paper rings with allergen	01 x 10
EP F111-0110	22312	<b>F111: Peppermint tea</b>	paper rings with allergen	01 x 10
EP F113-0110	22313	<b>F113: Black tea</b>	paper rings with allergen	01 x 10
EP F122-0110	22742	<b>F122: Olive, black</b>	paper rings with allergen	01 x 10
EP F126-0110	22611	<b>F126: Peppermint</b>	paper rings with allergen	01 x 10
EP F128-0110	22743	<b>F128: Poppy seed</b>	paper rings with allergen	01 x 10
EP F130-0110	22658	<b>F130: Turkey</b>	paper rings with allergen	01 x 10
EP F132-0110	22744	<b>F132: Bean, green</b>	paper rings with allergen	01 x 10
EP F142-0110	22745	<b>F142: Comfrey</b>	paper rings with allergen	01 x 10
EP F144-0110	22746	<b>F144: Pistachio</b>	paper rings with allergen	01 x 10
EP F145-0110	22747	<b>F145: Fig</b>	paper rings with allergen	01 x 10
EP F147-0110	22748	<b>F147: Hibiscus tea</b>	paper rings with allergen	01 x 10
EP F148-0110	22612	<b>F148: Plum</b>	paper rings with allergen	01 x 10
EP F149-0110	22749	<b>F149: Papaya</b>	paper rings with allergen	01 x 10
EP F150-0110	22613	<b>F150: Edam cheese</b>	paper rings with allergen	01 x 10
EP F151-0110	22750	<b>F151: Courgette (zucchini)</b>	paper rings with allergen	01 x 10
EP F152-0110	22614	<b>F152: Apricot</b>	paper rings with allergen	01 x 10
EP F153-0110	22615	<b>F153: Wheat bran</b>	paper rings with allergen	01 x 10
EP F154-0110	22751	<b>F154: St. Johns bread</b>	paper rings with allergen	01 x 10
EP F156-0110	22616	<b>F156: Raspberry</b>	paper rings with allergen	01 x 10
EP F158-0110	22617	<b>F158: Cashew nut</b>	paper rings with allergen	01 x 10
EP F159-0110	22752	<b>F159: Chickpea</b>	paper rings with allergen	01 x 10
EP F160-0110	22753	<b>F160: Sardine</b>	paper rings with allergen	01 x 10
EP F162-0110	22618	<b>F162: Dover sole</b>	paper rings with allergen	01 x 10
EP F163-0110	22619	<b>F163: Pike</b>	paper rings with allergen	01 x 10
EP F164-0110	22754	<b>F164: Swordfish</b>	paper rings with allergen	01 x 10

**Allercoat 6™-System: Allergen Rings (Food)**

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP F165-0110	22659	<b>F165: Veal, raw</b>	paper rings with allergen	01 x 10
EP F166-0110	22755	<b>F166: Chinese cabbage</b>	paper rings with allergen	01 x 10
EP F167-0110	22756	<b>F167: Rabbit meat, raw</b>	paper rings with allergen	01 x 10
EP F168-0110	22660	<b>F168: Milk powder</b>	paper rings with allergen	01 x 10
EP F170-0110	22620	<b>F170: Nectarine</b>	paper rings with allergen	01 x 10
EP F171-0110	22621	<b>F171: Currant, red</b>	paper rings with allergen	01 x 10
EP F172-0110	22758	<b>F172: Artichoke</b>	paper rings with allergen	01 x 10
EP F173-0110	22759	<b>F173: Gurnard</b>	paper rings with allergen	01 x 10
EP F175-0110	22622	<b>F175: Blackberry</b>	paper rings with allergen	01 x 10
EP F176-0110	22760	<b>F176: Aubergine (eggplant)</b>	paper rings with allergen	01 x 10
EP F177-0110	22761	<b>F177: Oyster</b>	paper rings with allergen	01 x 10
EP F178-0110	22623	<b>F178: Haddock</b>	paper rings with allergen	01 x 10
EP F180-0110	22763	<b>F180: Carp</b>	paper rings with allergen	01 x 10
EP F181-0110	22764	<b>F181: Shark</b>	paper rings with allergen	01 x 10
EP F182-0110	22765	<b>F182: Cranberry</b>	paper rings with allergen	01 x 10
EP F185-0110	22766	<b>F185: Cabbage, red</b>	paper rings with allergen	01 x 10
EP F187-0110	22661	<b>F187: Sheep's cheese</b>	paper rings with allergen	01 x 10
EP F188-0110	22624	<b>F188: Soy flour</b>	paper rings with allergen	01 x 10
EP F189-0110	22625	<b>F189: Rye bran</b>	paper rings with allergen	01 x 10
EP F190-0110	22626	<b>F190: Barley bran</b>	paper rings with allergen	01 x 10
EP F191-0110	22767	<b>F191: Pumpkin</b>	paper rings with allergen	01 x 10
EP F192-0110	22768	<b>F192: Quail</b>	paper rings with allergen	01 x 10
EP F193-0110	22769	<b>F193: Speit, unripe</b>	paper rings with allergen	01 x 10
EP F194-0110	22770	<b>F194: Pepperoni</b>	paper rings with allergen	01 x 10
EP F195-0110	22627	<b>F195: Soy milk</b>	paper rings with allergen	01 x 10
EP F196-0110	22662	<b>F196: Halibut</b>	paper rings with allergen	01 x 10
EP F197-0110	22628	<b>F197: Pine kernels/ nuts</b>	paper rings with allergen	01 x 10
EP F198-0110	22629	<b>F198: Gouda cheese</b>	paper rings with allergen	01 x 10
EP F199-0110	22771	<b>F199: Raisin</b>	paper rings with allergen	01 x 10
EP F201-0110	22772	<b>F201: Chanterelle, yellow</b>	paper rings with allergen	01 x 10
EP F203-0110	22773	<b>F203: Pomegranate</b>	paper rings with allergen	01 x 10
EP F205-0110	22630	<b>F205: Goat's cheese</b>	paper rings with allergen	01 x 10
EP F206-0110	22320	<b>F206: Mackerel</b>	paper rings with allergen	01 x 10
EP F209-0110	22775	<b>F209: Sage tea</b>	paper rings with allergen	01 x 10
EP F210-0110	22631	<b>F210: Wheat</b>	paper rings with allergen	01 x 10
EP F212-0110	22322	<b>F212: Mushroom</b>	paper rings with allergen	01 x 10
EP F213-0110	22776	<b>F213: Rhubarb</b>	paper rings with allergen	01 x 10
EP F214-0110	22777	<b>F214: Radish</b>	paper rings with allergen	01 x 10
EP F215-0110	22778	<b>F215: Starch (Maize)</b>	paper rings with allergen	01 x 10
EP F217-0110	22633	<b>F217: Soy protein</b>	paper rings with allergen	01 x 10

Allercoat 6™-System: Allergen Rings (Food)				
Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP F218-0110	22634	<b>F218: Lecithin</b>	paper rings with allergen	01 x 10
EP F219-0110	22635	<b>F219: Goat milk</b>	paper rings with allergen	01 x 10
EP F220-0110	22779	<b>F220: Anchovy</b>	paper rings with allergen	01 x 10
EP F221-0110	22819	<b>F221: Bamboo shoots</b>	paper rings with allergen	01 x 10
EP F222-0110	22780	<b>F222: Pumpkin seeds</b>	paper rings with allergen	01 x 10
EP F224-0110	22781	<b>F224: Beetroot</b>	paper rings with allergen	01 x 10
EP F226-0110	22782	<b>F226: Duck</b>	paper rings with allergen	01 x 10
EP F227-0110	22783	<b>F227: Venison</b>	paper rings with allergen	01 x 10
EP F228-0110	22784	<b>F228: Wild boar</b>	paper rings with allergen	01 x 10
EP F229-0110	22785	<b>F229: Blueberry</b>	paper rings with allergen	01 x 10
EP F230-0110	22820	<b>F230: Caviar</b>	paper rings with allergen	01 x 10
EP F231-0110	22325	<b>F231: Milk, cooked</b>	paper rings with allergen	01 x 10
EP F232-0110	22787	<b>F232: Monkfish</b>	paper rings with allergen	01 x 10
EP F233-0110	22788	<b>F233: Green kale</b>	paper rings with allergen	01 x 10
EP F234-0110	22326	<b>F234: Vanilla pod</b>	paper rings with allergen	01 x 10
EP F235-0110	22327	<b>F235: Lentil</b>	paper rings with allergen	01 x 10
EP F238-0110	22328	<b>F238: alpha-Lactalbumin, cooked</b>	paper rings with allergen	01 x 10
EP F239-0110	22329	<b>F239: beta-Lactoglobulin, cooked</b>	paper rings with allergen	01 x 10
EP F240-0110	22331	<b>F240: Casein, cooked</b>	paper rings with allergen	01 x 10
EP F244-0110	22332	<b>F244: Cucumber</b>	paper rings with allergen	01 x 10
EP F247-0110	22299	<b>F247: Honey</b>	paper rings with allergen	01 x 10
EP F254-0110	22333	<b>F254: Plaice</b>	paper rings with allergen	01 x 10
EP F258-0110	22334	<b>F258: Squid (Inkfish)</b>	paper rings with allergen	01 x 10
EP F260-0110	22335	<b>F260: Broccoli</b>	paper rings with allergen	01 x 10
EP F261-0110	22336	<b>F261: Asparagus</b>	paper rings with allergen	01 x 10
EP F265-0110	22337	<b>F265: Caraway</b>	paper rings with allergen	01 x 10
EP F267-0110	22338	<b>F267: Cardamom</b>	paper rings with allergen	01 x 10
EP F268-0110	22339	<b>F268: Clove</b>	paper rings with allergen	01 x 10
EP F269-0110	22340	<b>F269: Basil</b>	paper rings with allergen	01 x 10
EP F270-0110	22341	<b>F270: Ginger</b>	paper rings with allergen	01 x 10
EP F271-0110	22342	<b>F271: Anise</b>	paper rings with allergen	01 x 10
EP F272-0110	22343	<b>F272: Tarragon</b>	paper rings with allergen	01 x 10
EP F273-0110	22344	<b>F273: Thyme</b>	paper rings with allergen	01 x 10
EP F274-0110	22345	<b>F274: Marjoram</b>	paper rings with allergen	01 x 10
EP F275-0110	22346	<b>F275: Lovage</b>	paper rings with allergen	01 x 10
EP F276-0110	22347	<b>F276: Fennel</b>	paper rings with allergen	01 x 10
EP F281-0110	22348	<b>F281: Curry</b>	paper rings with allergen	01 x 10
EP F284-0110	22349	<b>F284: Turkey</b>	paper rings with allergen	01 x 10
EP F290-0125	22092	<b>F290: Milk</b>	paper rings with allergen	01 x 25
EP F291-0110	22351	<b>F291: Cauliflower</b>	paper rings with allergen	01 x 10

**Allercoat 6™-System: Allergen Rings (Food)**

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP F292-0110	22638	F292: Maize (Corn)	paper rings with allergen	01 x 10
EP F293-0110	22791	F293: Iceberg lettuce	paper rings with allergen	01 x 10
EP F294-0110	22352	F294: Grenadilla (Passion fruit)	paper rings with allergen	01 x 10
EP F296-0110	22794	F296: Gelatin, bovine	paper rings with allergen	01 x 10
EP F297-0110	22795	F297: Potato flour	paper rings with allergen	01 x 10
EP F298-0110	22796	F298: Coconut flakes	paper rings with allergen	01 x 10
EP F299-0110	22797	F299: Parmesan cheese	paper rings with allergen	01 x 10
EP F300-0110	22798	F300: Horse	paper rings with allergen	01 x 10
EP F301-0110	22799	F301: Sorrel	paper rings with allergen	01 x 10
EP F302-0110	22800	F302: Sauerkraut	paper rings with allergen	01 x 10
EP F303-0110	22801	F303: Snail	paper rings with allergen	01 x 10
EP F304-0110	22802	F304: Savoy cabbage	paper rings with allergen	01 x 10
EP F306-0110	22639	F306: Guar flour	paper rings with allergen	01 x 10
EP F307-0110	22640	F307: Bay leaf	paper rings with allergen	01 x 10
EP F308-0110	22641	F308: Cinnamon	paper rings with allergen	01 x 10
EP F309-0110	22642	F309: Coriander	paper rings with allergen	01 x 10
EP F310-0110	22643	F310: Oregano	paper rings with allergen	01 x 10
EP F311-0110	22644	F311: Peppercorn, green	paper rings with allergen	01 x 10
EP F312-0110	22645	F312: Pepper, white	paper rings with allergen	01 x 10
EP F313-0110	22803	F313: Chilli	paper rings with allergen	01 x 10
EP F314-0110	22804	F314: Savory	paper rings with allergen	01 x 10
EP F315-0110	22805	F315: Dill	paper rings with allergen	01 x 10
EP F316-0110	22806	F316: Chervil	paper rings with allergen	01 x 10
EP F317-0110	22807	F317: Balm	paper rings with allergen	01 x 10
EP F318-0110	22808	F318: Rosemary	paper rings with allergen	01 x 10
EP F319-0110	22809	F319: Sage	paper rings with allergen	01 x 10
EP F320-0110	22810	F320: Chive	paper rings with allergen	01 x 10
EP F321-0110	22811	F321: Juniper berry	paper rings with allergen	01 x 10
EP F322-0110	22656	F322: Ovalbumin	paper rings with allergen	01 x 10
EP F323-0110	22657	F323: Camembert	paper rings with allergen	01 x 10
EP F324-0110	22740	F324: Clam	paper rings with allergen	01 x 10
EP F325-0110	22774	F325: Bean, red	paper rings with allergen	01 x 10
EP F326-0110	22632	F326: Currant, black	paper rings with allergen	01 x 10
EP F327-0110	22637	F327: Wine vinegar	paper rings with allergen	01 x 10
EP F328-0110	22790	F328: Gooseberry	paper rings with allergen	01 x 10
EP F329-0110	22815	F329: Chamomile tea	paper rings with allergen	01 x 10
EP F330-0110	22816	F330: Linseed	paper rings with allergen	01 x 10
EP F331-0110	22741	F331: Brussels sprouts	paper rings with allergen	01 x 10
EP F332-0110	22786	F332: Lychee chinensis	paper rings with allergen	01 x 10
EP F333-0110	22789	F333: Chicory	paper rings with allergen	01 x 10

Allercoat 6™-System: Allergen Rings (Food)				
Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP F334-0110	22792	<b>F334: Endive</b>	paper rings with allergen	01 x 10
EP F335-0110	22814	<b>F335: Roquefort</b>	paper rings with allergen	01 x 10
EP F336-0110	22886	<b>F336: Glutamate</b>	paper rings with allergen	01 x 10
EP FS11-0125	22152	<b>FS11: Nuts (F13, F17, F18, F20, F36)</b>	paper rings with allergen mixture	01 x 25
EP FS12-0125	22153	<b>FS12: Seafood (F3, F24, F37, F40, F41)</b>	paper rings with allergen mixture	01 x 25
EP FS13-0125	22154	<b>FS13: Mixed cereals (F4, F7, F8, F10, F11)</b>	paper rings with allergen mixture	01 x 25
EP FS14-0110	22981	<b>FS14: Fish (F3, F21, F206, F254)</b>	paper rings with allergen mixture	01 x 10
EP FS15-0125	22155	<b>FS15: Baby food (F1, F2, F4, F13, F14)</b>	paper rings with allergen mixture	01 x 25
EP FS16-0110	22951	<b>FS16: Food I (F15, F12, F235, F14, F13)</b>	paper rings with allergen mixture	01 x 10
EP FS17-0110	22952	<b>FS17: Food II (F4, F5, F6, F99)</b>	paper rings with allergen mixture	01 x 10
EP FS18-0110	22953	<b>FS18: Food III (F11, F8, F18, F101, F104)</b>	paper rings with allergen mixture	01 x 10
EP FS19-0110	22954	<b>FS19: Food IV (F49, F94, F95, F97)</b>	paper rings with allergen mixture	01 x 10
EP FS20-0110	22955	<b>FS20: Food V (safflower oil, U87, F100, F272)</b>	paper rings with allergen mixture	01 x 10
EP FS21-0125	22169	<b>FS21: Fruit I (F29, F33, F95, F49)</b>	paper rings with allergen mixture	01 x 25
EP FS22-0125	22170	<b>FS22: Fruit II (F32, F44, F72, F94)</b>	paper rings with allergen mixture	01 x 25
EP FS23-0125	22172	<b>FS23: Flours I (F7, F5, F4, F79)</b>	paper rings with allergen mixture	01 x 25
EP FS24-0125	22173	<b>FS24: Vegetables I (F12, F15, F31, F35)</b>	paper rings with allergen mixture	01 x 25
EP FS25-0125	22174	<b>FS25: Vegetables II (F25, F39, F46, F38)</b>	paper rings with allergen mixture	01 x 25
EP FS26-0110	22964	<b>FS26: Vegetables III (F14, F85, F108, F212)</b>	paper rings with allergen mixture	01 x 10
EP FS27-0110	22968	<b>FS27: Cheese (F70, F81, F82, F150)</b>	paper rings with allergen mixture	01 x 10
EP FS28-0110	22970	<b>FS28: Meat (F26, F27, F88)</b>	paper rings with allergen mixture	01 x 10
EP FS29-0110	22975	<b>FS29: Poultry (F57, F58, F83, F130)</b>	paper rings with allergen mixture	01 x 10
EP FS30-0110	22976	<b>FS30: Flours II (F4, F5, F6, F7, F8)</b>	paper rings with allergen mixture	01 x 10
EP FS31-0110	22979	<b>FS31: Food VI (F3, F8, F17, F23, F25)</b>	paper rings with allergen mixture	01 x 10
EP FS32-0110	22982	<b>FS32: Citrus fruits (F30, F32, F33, F34)</b>	paper rings with allergen mixture	01 x 10
EP FS55-0125	22178	<b>FS55: Multiscreen soy I (F14, F51, F188, F195, F217)</b>	paper rings with allergen mixture	01 x 25
EP FS70-0110	22956	<b>FS70: Spices I (F272, F274, F273, F275)</b>	paper rings with allergen mixture	01 x 10

**Allercoat 6™-System: Allergen Rings (Food)**

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP FS71-0125	22156	FS71: Spices II (F265, F268, F92, F267)	paper rings with allergen mixture	01 x 25
EP FS72-0110	22957	FS72: Spices III (F271, F270, F269, F276)	paper rings with allergen mixture	01 x 10
EP FS73-0110	22958	FS73: Spices IV (F275, F265, F276, F271)	paper rings with allergen mixture	01 x 10
EP FS74-0110	22959	FS74: Spices V (F274, F273, F269, F111)	paper rings with allergen mixture	01 x 10
EP FS75-0125	22175	FS75: Spices VI (F271, F281, F265, F47)	paper rings with allergen mixture	01 x 25
EP FS76-0125	22176	FS76: Spices VII (F89, F90, F91, F307)	paper rings with allergen mixture	01 x 25

Allercoat™ 6-System: Allergen Rings (Grasses)				
Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP G001-0110	22396	<b>G1: Sweet vernal grass</b>	paper rings with allergen	01 x 10
EP G002-0110	22397	<b>G2: Bermuda grass</b>	paper rings with allergen	01 x 10
EP G003-0110	22398	<b>G3: Cocksfoot</b>	paper rings with allergen	01 x 10
EP G004-0110	22399	<b>G4: Fescue, meadow</b>	paper rings with allergen	01 x 10
EP G005-0125	22076	<b>G5: Rye grass, perennial</b>	paper rings with allergen	01 x 25
EP G006-0125	22077	<b>G6: Timothy grass</b>	paper rings with allergen	01 x 25
EP G008-0110	22401	<b>G8: Kentucky blue (June) grass</b>	paper rings with allergen	01 x 10
EP G009-0110	22402	<b>G9: Bentgrass</b>	paper rings with allergen	01 x 10
EP G010-0110	22403	<b>G10: Johnson Grass</b>	paper rings with allergen	01 x 10
EP G011-0110	22406	<b>G11: Brome, smooth</b>	paper rings with allergen	01 x 10
EP G012-0125	22078	<b>G12: Rye, cultivated</b>	paper rings with allergen	01 x 25
EP G013-0110	22407	<b>G13: Velvet grass</b>	paper rings with allergen	01 x 10
EP G014-0110	22408	<b>G14: Oats, cultivated</b>	paper rings with allergen	01 x 10
EP G016-0110	22409	<b>G16: Meadow foxtail</b>	paper rings with allergen	01 x 10
EP G017-0110	22410	<b>G17: Bahia grass</b>	paper rings with allergen	01 x 10
EP G018-0110	22411	<b>G18: Barley, cultivated</b>	paper rings with allergen	01 x 10
EP G019-0110	22821	<b>G19: Dogtail grass, crested</b>	paper rings with allergen	01 x 10
EP G020-0110	22412	<b>G20: Maize (Corn)</b>	paper rings with allergen	01 x 10
EP G021-0110	22828	<b>G21: Twich (Couch) grass</b>	paper rings with allergen	01 x 10
EP G070-0110	22413	<b>G70: Elymus</b>	paper rings with allergen	01 x 10
EP G071-0110	22822	<b>G71: Oat grass</b>	paper rings with allergen	01 x 10
EP G072-0110	22823	<b>G72: Knot grass</b>	paper rings with allergen	01 x 10
EP G073-0110	22824	<b>G73: Elymus triticoides</b>	paper rings with allergen	01 x 10
EP G074-0110	22825	<b>G74: Phalaris arundinacea</b>	paper rings with allergen	01 x 10
EP G075-0110	22826	<b>G75: Hordeum leporinum</b>	paper rings with allergen	01 x 10
EP G076-0110	22827	<b>G76: Wild oat</b>	paper rings with allergen	01 x 10
EP G100-0110	22404	<b>G100: Russian thistle (Saltwort)</b>	paper rings with allergen	01 x 10
EP G101-0110	22405	<b>G101: Wheat</b>	paper rings with allergen	01 x 10
EP GIH1-0125	22161	<b>IHS1: Inhalation I (G6, W6, W9, W21, T3)</b>	paper rings with allergen mixture	01 x 25
EP GIH2-0125	22162	<b>IHS2: Inhalation II (D2, E1, E2, E3, M6)</b>	paper rings with allergen mixture	01 x 25
EP GIH3-0110	22960	<b>IHS3: Inhalation III (D1, D2, E1, E2, M5)</b>	paper rings with allergen mixture	01 x 10
EP GIH4-0125	22179	<b>IHS4: Inhalation IV (G2, G5, G17, W1, W9, W10)</b>	paper rings with allergen mixture	01 x 25
EP GS01-0125	22158	<b>GS1: Grass (G1, G2, G5, G6, G10)</b>	paper rings with allergen mixture	01 x 25
EP GS11-0125	22157	<b>GS11: Grass, early flowering (G3, G4, G5, G6, G8)</b>	paper rings with allergen mixture	01 x 25
EP GS14-0125	22159	<b>GS14: Grass, late flowering (G1, G5, G7, G12, G13)</b>	paper rings with allergen mixture	01 x 25

## Allercoat™ 6-System: Allergen Rings (Grasses)

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP GS20-0125	22177	GS20: Cereal pollen (G12, G14, G18, G20, G101)	paper rings with allergen mixture	01 x 25
EP GS22-0125	22180	GS22: Grass Mix Australia (G2, G5, G6, G8, G10, G17)	paper rings with allergen mixture	01 x 25

## Allercoat 6™-System: Allergen Rings (House Dust)

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP H002-0125	22080	H2: House dust Hollister-Stier	paper rings with allergen	01 x 25
EP HS12-0125	22160	HS12: Mites (H2, D1, D2, I6)	paper rings with allergen mixture	01 x 25

## Allercoat™ 6-System: Allergen Rings (Insects)

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP I001-0110	22415	I1: Honeybee, venom	paper rings with allergen	01 x 10
EP I002-0110	22416	I2: White faced hornet ( <i>Dolichovespula maculata</i> )	paper rings with allergen	01 x 10
EP I003-0110	22417	I3: Wasp (Yellow jacket), venom	paper rings with allergen	01 x 10
EP I004-0110	22418	I4: Paper wasp ( <i>Polistes apachus</i> )	paper rings with allergen	01 x 10
EP I005-0110	22419	I5: Yellow hornet ( <i>Dolichovespula arenaria</i> )	paper rings with allergen	01 x 10
EP I006-0110	22420	I6: German cockroach	paper rings with allergen	01 x 10
EP I007-0110	22421	I7: Bumble bee, venom	paper rings with allergen	01 x 10
EP I010-0110	22829	I10: Paper wasp, venom	paper rings with allergen	01 x 10
EP I011-0110	22830	I11: Phospholipase	paper rings with allergen	01 x 10
EP I012-0110	22422	I12: Black ant	paper rings with allergen	01 x 10
EP I014-0110	22423	I14: Horsefly	paper rings with allergen	01 x 10
EP I015-0110	22831	I15: House fly	paper rings with allergen	01 x 10
EP I070-0110	22424	I70: Fire ant	paper rings with allergen	01 x 10
EP I071-0110	22425	I71: Mosquito (Gnat)	paper rings with allergen	01 x 10
EP I073-0110	22427	I73: Red gnat larvae	paper rings with allergen	01 x 10
EP I074-0110	22832	I74: Waterflea	paper rings with allergen	01 x 10
EP I075-0110	22428	I75: Hornet, venom	paper rings with allergen	01 x 10
EP I076-0110	22569	I76: <i>Trogoderma angustum</i>	paper rings with allergen	01 x 10
EP I077-0110	22429	I77: Rice flour beetle	paper rings with allergen	01 x 10

## Allercoat™ 6-System: Allergen Rings (Moulds)

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP M001-0125	22081	M1: <i>Penicillium notatum</i>	paper rings with allergen	01 x 25
EP M002-0125	22082	M2: <i>Cladosporium herbarum</i>	paper rings with allergen	01 x 25
EP M003-0125	22083	M3: <i>Aspergillus fumigatus</i>	paper rings with allergen	01 x 25
EP M004-0110	22430	M4: <i>Mucor racemosus</i>	paper rings with allergen	01 x 10
EP M005-0125	22084	M5: <i>Candida albicans</i>	paper rings with allergen	01 x 25
EP M006-0125	22085	M6: <i>Alternaria tenuis</i>	paper rings with allergen	01 x 25
EP M007-0110	22431	M7: <i>Botrytis cinerea</i>	paper rings with allergen	01 x 10
EP M008-0110	22432	M8: <i>Helminthosporium halodes</i>	paper rings with allergen	01 x 10
EP M009-0110	22433	M9: <i>Fusarium culmorum</i>	paper rings with allergen	01 x 10
EP M010-0110	22434	M10: <i>Stemphylium botryosum</i>	paper rings with allergen	01 x 10
EP M011-0110	22435	M11: <i>Rhizopus nigricans</i>	paper rings with allergen	01 x 10
EP M012-0110	22436	M12: <i>Aureobasidium pullulans</i>	paper rings with allergen	01 x 10
EP M013-0110	22437	M13: <i>Phoma betae</i>	paper rings with allergen	01 x 10
EP M014-0110	22438	M14: <i>Epicoccum purpurascens</i>	paper rings with allergen	01 x 10
EP M015-0110	22439	M15: <i>Trichoderma viride</i>	paper rings with allergen	01 x 10
EP M016-0110	22440	M16: <i>Curvularia lunata</i>	paper rings with allergen	01 x 10
EP M017-0110	22441	M17: <i>Aspergillus amstelodami</i>	paper rings with allergen	01 x 10
EP M018-0110	22442	M18: <i>Aspergillus terreus</i>	paper rings with allergen	01 x 10
EP M019-0110	22443	M19: <i>Aspergillus versicolor</i>	paper rings with allergen	01 x 10
EP M020-0110	22444	M20: <i>Neurospora sitophila</i>	paper rings with allergen	01 x 10
EP M021-0110	22446	M21: <i>Paecilomyces spec.</i>	paper rings with allergen	01 x 10
EP M023-0110	22448	M23: <i>Serpula lacrymans</i>	paper rings with allergen	01 x 10
EP M025-0110	22450	M25: <i>Aspergillus flavus</i>	paper rings with allergen	01 x 10
EP M030-0110	22836	M30: <i>Penicillium roqueforti</i>	paper rings with allergen	01 x 10
EP M031-0110	22837	M31: <i>Penicillium viridicatum</i>	paper rings with allergen	01 x 10
EP M032-0110	22668	M32: <i>Cladosporium cladosporioides</i>	paper rings with allergen	01 x 10
EP M037-0110	22838	M37: <i>Trichophyton mentagrophytes</i>	paper rings with allergen	01 x 10
EP M038-0110	22839	M38: <i>Trichophyton verrucosum</i>	paper rings with allergen	01 x 10
EP M039-0110	22840	M39: <i>Trichophyton rubrum</i>	paper rings with allergen	01 x 10
EP M040-0110	22663	M40: <i>Fusarium moniliforme</i>	paper rings with allergen	01 x 10
EP M041-0110	22664	M41: <i>Mucor mucedo</i>	paper rings with allergen	01 x 10
EP M046-0110	22666	M46: <i>Mucor spinosus</i>	paper rings with allergen	01 x 10
EP M047-0110	22841	M47: <i>Cephalosporium acremonium</i>	paper rings with allergen	01 x 10
EP M048-0110	22842	M48: <i>Aspergillus oryzae</i>	paper rings with allergen	01 x 10
EP M049-0110	22843	M49: <i>Aspergillus nidulans</i>	paper rings with allergen	01 x 10
EP M056-0110	22847	M56: <i>Microsporum canis</i>	paper rings with allergen	01 x 10
EP M058-0110	22833	M58: <i>Cladosporium fulvum</i>	paper rings with allergen	01 x 10
EP M059-0110	22834	M59: <i>Penicillium commune</i>	paper rings with allergen	01 x 10
EP M060-0110	22885	M60: <i>Ulocladium chartarum</i>	paper rings with allergen	01 x 10
EP M061-0110	22835	M61: <i>Penicillium expansum</i>	paper rings with allergen	01 x 10

## Allercoat™ 6-System: Allergen Rings (Moulds)

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP M062-0110	22670	<b>M62: Wine yeast</b>	paper rings with allergen	01 x 10
EP M063-0110	22844	<b>M63: Chaetomium globosum</b>	paper rings with allergen	01 x 10
EP M064-0110	22845	<b>M64: Penicillium chrysogenum</b>	paper rings with allergen	01 x 10
EP M207-0110	22445	<b>M207: Aspergillus niger</b>	paper rings with allergen	01 x 10
EP MS01-0125	22163	<b>MS1: Moulds (M1, M2, M3, M6)</b>	paper rings with allergen mixture	01 x 25
EP MS02-0110	22966	<b>MS2: Aspergilli 1 (M3, M21, M17, M40)</b>	paper rings with allergen mixture	01 x 10
EP MS03-0110	22967	<b>MS3: Aspergilli 2 (M19, M28, M207, M18)</b>	paper rings with allergen mixture	01 x 10
EP MS04-0110	22969	<b>MS4: Penicillia (M1, M22, M28, M30)</b>	paper rings with allergen mixture	01 x 10

## Allercoat™ 6-System: Allergen Rings (Further Allergens)

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP O070-0110	22454	<b>O70: Seminal fluid</b>	paper rings with allergen	01 x 10
EP O071-0110	22850	<b>O71: Seminal plasma</b>	paper rings with allergen	01 x 10
EP O072-0110	22851	<b>O72: Sperm sediment</b>	paper rings with allergen	01 x 10
EP O201-0110	22453	<b>O201: Tobacco dust</b>	paper rings with allergen	01 x 10

Allercoat™ 6-System: Allergen Rings (Trees)				
Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP T001-0110	22458	T1: Maple	paper rings with allergen	01 x 10
EP T002-0125	22086	T2: Alder	paper rings with allergen	01 x 25
EP T003-0125	22087	T3: Birch	paper rings with allergen	01 x 25
EP T004-0125	22088	T4: Hazel	paper rings with allergen	01 x 25
EP T005-0125	22089	T5: Beech, European	paper rings with allergen	01 x 25
EP T006-0110	22459	T6: Cedar, Japanese	paper rings with allergen	01 x 10
EP T007-0110	22460	T7: Oak	paper rings with allergen	01 x 10
EP T008-0110	22461	T8: Elm	paper rings with allergen	01 x 10
EP T009-0110	22462	T9: Olive	paper rings with allergen	01 x 10
EP T010-0110	22463	T10: Walnut	paper rings with allergen	01 x 10
EP T011-0110	22470	T11: Planetree	paper rings with allergen	01 x 10
EP T012-0110	22471	T12: Willow	paper rings with allergen	01 x 10
EP T013-0110	22853	T13: Mock, orange	paper rings with allergen	01 x 10
EP T014-0110	22472	T14: Poplar	paper rings with allergen	01 x 10
EP T015-0110	22473	T15: Ash	paper rings with allergen	01 x 10
EP T016-0110	22474	T16: Pine	paper rings with allergen	01 x 10
EP T017-0110	22475	T17: Chestnut, European	paper rings with allergen	01 x 10
EP T018-0110	22476	T18: Eucalyptus	paper rings with allergen	01 x 10
EP T019-0110	22477	T19: Mimosa	paper rings with allergen	01 x 10
EP T020-0110	22478	T20: Privet, common	paper rings with allergen	01 x 10
EP T021-0110	22298	T21: Lilac, common	paper rings with allergen	01 x 10
EP T022-0110	22854	T22: Hawthorn	paper rings with allergen	01 x 10
EP T023-0110	22479	T23: Cypress	paper rings with allergen	01 x 10
EP T024-0110	22480	T24: Lime tree	paper rings with allergen	01 x 10
EP T025-0110	22481	T25: Elder, European	paper rings with allergen	01 x 10
EP T026-0110	22672	T26: Golden chain	paper rings with allergen	01 x 10
EP T027-0110	22852	T27: Cedar, mountain	paper rings with allergen	01 x 10
EP T028-0110	22673	T28: Robinia	paper rings with allergen	01 x 10
EP T029-0110	22855	T29: Morello cherry	paper rings with allergen	01 x 10
EP T031-0110	22856	T31: Melaleuca	paper rings with allergen	01 x 10
EP T032-0110	22857	T32: Orange	paper rings with allergen	01 x 10
EP T033-0110	22858	T33: Poplar, lombardy	paper rings with allergen	01 x 10
EP T037-0110	22859	T37: Yew tree	paper rings with allergen	01 x 10
EP T070-0110	22482	T70: Mulberry, white	paper rings with allergen	01 x 10
EP T100-0110	22464	T100: Pecan	paper rings with allergen	01 x 10
EP T101-0110	22465	T101: Walnut, black	paper rings with allergen	01 x 10
EP T102-0110	22466	T102: Mesquite	paper rings with allergen	01 x 10
EP T103-0110	22467	T103: Live Oak	paper rings with allergen	01 x 10
EP T104-0110	22468	T104: Australian pine	paper rings with allergen	01 x 10
EP T107-0110	22469	T107: Cajeput tree	paper rings with allergen	01 x 10

**Allercoat™ 6-System: Allergen Rings (Trees)**

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP TS01-0110	22961	TS1: Trees I (T3, T6, T8, T102, T103)	paper rings with allergen mixture	01 x 10
EP TS02-0110	22962	TS2: Trees II (T1, T3, T5, T11, T14, T103)	paper rings with allergen mixture	01 x 10
EP TS15-0125	22164	TS15: Trees, early flowering (T2, T4, T8, T12, T14)	paper rings with allergen mixture	01 x 25
EP TS16-0125	22165	TS16: Trees, late flowering (T1, T3, T5, T7, T10)	paper rings with allergen mixture	01 x 25
EP TS17-0110	22971	TS17: Soft woods (WD02, WD14, WD06, WD25)	paper rings with allergen mixture	01 x 10
EP TS18-0110	22972	TS18: Hard woods (WD20, WD17, WD12, WD9)	paper rings with allergen mixture	01 x 10
EP TS19-0125	22181	TS19: Tree Mix Australia (T104, T12, T18, T19, T107)	paper rings with allergen mixture	01 x 25

**Allercoat™ 6-System: Allergen Rings (Environmental Allergens)**

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP U070-0110	22483	U70: Green coffee bean	paper rings with allergen	01 x 10
EP U072-0110	22485	U72: Ispaghula ( <i>Sem. plant. psyll.</i> )	paper rings with allergen	01 x 10
EP U073-0110	22486	U73: Wild silk	paper rings with allergen	01 x 10
EP U074-0110	22487	U74: Raw silk, fabric	paper rings with allergen	01 x 10
EP U075-0110	22488	U75: TDI Isocyanate - HSA	paper rings with allergen	01 x 10
EP U076-0110	22489	U76: MDI Isocyanate - HSA	paper rings with allergen	01 x 10
EP U077-0110	22490	U77: HDI Isocyanate - HSA	paper rings with allergen	01 x 10
EP U078-0110	22491	U78: Ethylene oxide - HSA	paper rings with allergen	01 x 10
EP U079-0110	22492	U79: Phthalic anhydride - HSA	paper rings with allergen	01 x 10
EP U080-0110	22493	U80: Hops	paper rings with allergen	01 x 10
EP U081-0110	22494	U81: Cotton, yarn and fabric	paper rings with allergen	01 x 10
EP U082-0110	22495	U82: Sheep's wool	paper rings with allergen	01 x 10
EP U083-0110	22496	U83: Formaldehyde - HSA	paper rings with allergen	01 x 10
EP U084-0110	22497	U84: Ficus benjamina	paper rings with allergen	01 x 10
EP U085-0110	22498	U85: Latex	paper rings with allergen	01 x 10
EP U086-0110	22499	U86: Trimellitic acid anhydride (TMA - HSA)	paper rings with allergen	01 x 10
EP U087-0110	22500	U87: Sunflower seeds	paper rings with allergen	01 x 10
EP U088-0110	22501	U88: Chloramine T - HSA	paper rings with allergen	01 x 10
EP U089-0110	22502	U89: Alpha-amylase	paper rings with allergen	01 x 10
EP U091-0110	22504	U91: PCP - HSA	paper rings with allergen	01 x 10
EP U092-0110	22505	U92: Pyrethroide - HSA	paper rings with allergen	01 x 10
EP U093-0110	22538	U93: Maleic anhydride - HSA	paper rings with allergen	01 x 10

## Allercoat™ 6-System: Allergen Rings (Environmental Allergens)

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP U094-0110	22506	<b>U94: Nickel sulphate - HSA</b>	paper rings with allergen	01 x 10
EP U101-0110	22679	<b>U101: Benzoic acid - HSA</b>	paper rings with allergen	01 x 10
EP U102-0110	22680	<b>U102: Parabene propylester - HSA</b>	paper rings with allergen	01 x 10
EP U103-0110	22681	<b>U103: Sorbic acid - HSA</b>	paper rings with allergen	01 x 10
EP U106-0110	22709	<b>U106: Benzoic acid ester - HSA</b>	paper rings with allergen	01 x 10
EP U107-0110	22683	<b>U107: Butyylester - HSA</b>	paper rings with allergen	01 x 10
EP U108-0110	22684	<b>U108: Mercury - HSA</b>	paper rings with allergen	01 x 10
EP U109-0110	22685	<b>U109: Chromium - HSA</b>	paper rings with allergen	01 x 10
EP U110-0110	22686	<b>U110: Copper - HSA</b>	paper rings with allergen	01 x 10
EP U111-0110	22687	<b>U111: Cobalt - HSA</b>	paper rings with allergen	01 x 10
EP U112-0110	22688	<b>U112: Silver - HSA</b>	paper rings with allergen	01 x 10
EP U113-0110	22689	<b>U113: Gold - HSA</b>	paper rings with allergen	01 x 10
EP U115-0110	22691	<b>U115: Tin - HSA</b>	paper rings with allergen	01 x 10
EP U116-0110	22692	<b>U116: Zinc - HSA</b>	paper rings with allergen	01 x 10
EP U117-0110	22693	<b>U117: Formic acid - HSA</b>	paper rings with allergen	01 x 10
EP U118-0110	22889	<b>U118: Chinoline yellow - HSA</b>	paper rings with allergen	01 x 10
EP U119-0110	22890	<b>U119: Yellow orange - HSA</b>	paper rings with allergen	01 x 10
EP U121-0110	22697	<b>U121: Parabenemethylester - HSA</b>	paper rings with allergen	01 x 10
EP U122-0110	22698	<b>U122: Tartrazine - HSA</b>	paper rings with allergen	01 x 10
EP U124-0110	22700	<b>U124: Palladium - HSA</b>	paper rings with allergen	01 x 10
EP U125-0110	22701	<b>U125: Lanolin</b>	paper rings with allergen	01 x 10
EP U127-0110	22703	<b>U127: Rayon</b>	paper rings with allergen	01 x 10
EP U128-0110	22704	<b>U128: Linen</b>	paper rings with allergen	01 x 10
EP U129-0110	22706	<b>U129: Jute</b>	paper rings with allergen	01 x 10
EP U131-0110	22864	<b>U131: Acrilan</b>	paper rings with allergen	01 x 10
EP U132-0110	22865	<b>U132: Threshing dust</b>	paper rings with allergen	01 x 10
EP U133-0110	22866	<b>U133: Flax</b>	paper rings with allergen	01 x 10
EP U134-0110	22867	<b>U134: Ceiba pentandra</b>	paper rings with allergen	01 x 10
EP U135-0110	22868	<b>U135: Nylon</b>	paper rings with allergen	01 x 10
EP U136-0110	22869	<b>U136: Silk</b>	paper rings with allergen	01 x 10
EP U137-0110	22870	<b>U137: Terylene</b>	paper rings with allergen	01 x 10
EP U138-0110	22871	<b>U138: Wheat threshing dust</b>	paper rings with allergen	01 x 10
EP U140-0110	22708	<b>U140: Sheep's wool, raw</b>	paper rings with allergen	01 x 10
EP U141-0110	22705	<b>U141: Hemicellulase</b>	paper rings with allergen	01 x 10
EP U142-0110	22544	<b>U142: Alcalase</b>	paper rings with allergen	01 x 10
EP U145-0110	22817	<b>U145: Henna, red</b>	paper rings with allergen	01 x 10
EP U146-0110	22818	<b>U146: Collagen</b>	paper rings with allergen	01 x 10
EP US11-0125	22166	<b>US11: Environmental screen I (E1, E2, D1, D2, H2)</b>	paper rings with allergen mixture	01 x 25

**Allercoat™ 6-System: Allergen Rings (Environmental Allergens)**

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP US12-0110	22974	<b>US12: Environmental screen II (G12, W6, T3, E1, E5, D1, M2)</b>	paper rings with allergen mixture	01 x 10
EP US13-0110	22977	<b>US13: Synthetic fibres (U131, U127, U135, U137)</b>	paper rings with allergen mixture	01 x 10
EP US14-0110	22978	<b>US14: Natural fibres (U82, U81, U129, U74)</b>	paper rings with allergen mixture	01 x 10
EP US15-0110	22980	<b>US15: Chemicals I (U75, U76, U77, U79)</b>	paper rings with allergen mixture	01 x 10

**Allercoat™ 6-System: Allergen Rings (Herbal and Flower Pollen)**

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP W001-0110	22507	<b>W1: Ragweed, short</b>	paper rings with allergen	01 x 10
EP W002-0110	22508	<b>W2: Ragweed, western</b>	paper rings with allergen	01 x 10
EP W003-0110	22509	<b>W3: Ragweed, tall or giant</b>	paper rings with allergen	01 x 10
EP W004-0110	22510	<b>W4: Ragweed, false</b>	paper rings with allergen	01 x 10
EP W005-0110	22511	<b>W5: Wormwood</b>	paper rings with allergen	01 x 10
EP W006-0125	22090	<b>W6: Mugwort</b>	paper rings with allergen	01 x 25
EP W007-0110	22512	<b>W7: Ox-eye daisy (Marguerite)</b>	paper rings with allergen	01 x 10
EP W008-0110	22513	<b>W8: Dandelion</b>	paper rings with allergen	01 x 10
EP W009-0125	22091	<b>W9: Plantain, English</b>	paper rings with allergen	01 x 25
EP W010-0110	22514	<b>W10: Lamb's quarter</b>	paper rings with allergen	01 x 10
EP W011-0110	22519	<b>W11: Russian thistle (Saltwort)</b>	paper rings with allergen	01 x 10
EP W012-0110	22520	<b>W12: Goldenrod</b>	paper rings with allergen	01 x 10
EP W013-0110	22521	<b>W13: Cocklebur</b>	paper rings with allergen	01 x 10
EP W014-0110	22522	<b>W14: Amaranth</b>	paper rings with allergen	01 x 10
EP W015-0110	22523	<b>W15: Lens-scale</b>	paper rings with allergen	01 x 10
EP W016-0110	22524	<b>W16: Rough marsh/ Elder</b>	paper rings with allergen	01 x 10
EP W017-0110	22525	<b>W17: Firebush</b>	paper rings with allergen	01 x 10
EP W019-0110	22526	<b>W19: Wall pellitory (Parietaria officinalis)</b>	paper rings with allergen	01 x 10
EP W021-0110	22527	<b>W21: Wall pellitory, erect (Parietaria judaica)</b>	paper rings with allergen	01 x 10
EP W022-0110	22710	<b>W22: Chrysanthemum</b>	paper rings with allergen	01 x 10
EP W023-0110	22712	<b>W23: Dahlia</b>	paper rings with allergen	01 x 10
EP W024-0110	22874	<b>W24: Kochia scoparia</b>	paper rings with allergen	01 x 10
EP W025-0110	22528	<b>W25: Geranium</b>	paper rings with allergen	01 x 10

Allercoat™ 6-System: Allergen Rings (Herbal and Flower Pollen)				
Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP W026-0110	22713	<b>W26: Daffodil</b>	paper rings with allergen	01 x 10
EP W027-0110	22714	<b>W27: Carnation, whole flower</b>	paper rings with allergen	01 x 10
EP W028-0110	22529	<b>W28: Rose</b>	paper rings with allergen	01 x 10
EP W029-0110	22875	<b>W29: Sunflower</b>	paper rings with allergen	01 x 10
EP W030-0110	22530	<b>W30: Tulip</b>	paper rings with allergen	01 x 10
EP W031-0110	22531	<b>W31: Rape</b>	paper rings with allergen	01 x 10
EP W032-0110	22532	<b>W32: Chamomile, wild</b>	paper rings with allergen	01 x 10
EP W033-0110	22533	<b>W33: Primrose</b>	paper rings with allergen	01 x 10
EP W034-0110	22884	<b>W34: Clover, red</b>	paper rings with allergen	01 x 10
EP W038-0110	22877	<b>W38: Marchelder, rough</b>	paper rings with allergen	01 x 10
EP W039-0110	22715	<b>W39: Lupine</b>	paper rings with allergen	01 x 10
EP W040-0110	22716	<b>W40: Hyacinth</b>	paper rings with allergen	01 x 10
EP W041-0110	22878	<b>W41: Alfalfa (lucerne)</b>	paper rings with allergen	01 x 10
EP W042-0110	22879	<b>W42: Lady's mantle</b>	paper rings with allergen	01 x 10
EP W044-0110	22717	<b>W44: Lily</b>	paper rings with allergen	01 x 10
EP W045-0110	22880	<b>W45: Euphorbia</b>	paper rings with allergen	01 x 10
EP W046-0110	22718	<b>W46: Azalea</b>	paper rings with allergen	01 x 10
EP W047-0110	22881	<b>W47: Hibiscus</b>	paper rings with allergen	01 x 10
EP W048-0110	22882	<b>W48: Yarrow</b>	paper rings with allergen	01 x 10
EP W049-0110	22719	<b>W49: Begonia, flowers and leaves</b>	paper rings with allergen	01 x 10
EP W050-0110	22720	<b>W50: Forsythia</b>	paper rings with allergen	01 x 10
EP W051-0110	22711	<b>W51: Aster</b>	paper rings with allergen	01 x 10
EP W052-0110	22873	<b>W52: Willow herb</b>	paper rings with allergen	01 x 10
EP W053-0110	22876	<b>W53: Heather</b>	paper rings with allergen	01 x 10
EP W054-0110	22883	<b>W54: Coltsfoot</b>	paper rings with allergen	01 x 10
EP W100-0110	22515	<b>W100: Sorrel, red</b>	paper rings with allergen	01 x 10
EP W101-0110	22516	<b>W101: Spiny pigweed</b>	paper rings with allergen	01 x 10
EP W102-0110	22517	<b>W102: Careless weed</b>	paper rings with allergen	01 x 10
EP W103-0110	22518	<b>W103: Nettle</b>	paper rings with allergen	01 x 10
EP WD01-0110	22534	<b>WD1: Sycamore maple</b>	paper rings with allergen	01 x 10
EP WD02-0110	22535	<b>WD2: Beech, European</b>	paper rings with allergen	01 x 10
EP WD03-0110	22536	<b>WD3: Oak</b>	paper rings with allergen	01 x 10
EP WD04-0110	22537	<b>WD4: Ash</b>	paper rings with allergen	01 x 10
EP WD05-0110	22675	<b>WD5: Spruce, European</b>	paper rings with allergen	01 x 10
EP WD06-0110	22539	<b>WD6: Pine</b>	paper rings with allergen	01 x 10
EP WD07-0110	22540	<b>WD7: Limba (Afara)</b>	paper rings with allergen	01 x 10
EP WD08-0110	22541	<b>WD8: Makore</b>	paper rings with allergen	01 x 10
EP WD09-0110	22542	<b>WD9: Mahogany, African</b>	paper rings with allergen	01 x 10
EP WD12-0110	22543	<b>WD12: Ramin</b>	paper rings with allergen	01 x 10
EP WD14-0110	22677	<b>WD14: Fir</b>	paper rings with allergen	01 x 10

**Allercoat™ 6-System: Allergen Rings (Herbal and Flower Pollen)**

Order No.	Previous Order No.	Antibodies of Class IgE against	Substrate	Format
EP WD15-0110	22545	<b>WD15: Teak</b>	paper rings with allergen	01 x 10
EP WD16-0110	22546	<b>WD16: Afzelia</b>	paper rings with allergen	01 x 10
EP WD18-0110	22547	<b>WD18: Cherry wood</b>	paper rings with allergen	01 x 10
EP WD19-0110	22548	<b>WD19: Meranti (Red Lauan)</b>	paper rings with allergen	01 x 10
EP WD21-0110	22676	<b>WD21: Acacia</b>	paper rings with allergen	01 x 10
EP WD22-0110	22678	<b>WD22: Sapeli</b>	paper rings with allergen	01 x 10
EP WD23-0110	22863	<b>WD23: Poplar wood</b>	paper rings with allergen	01 x 10
EP WD24-0110	22862	<b>WD24: Lime wood</b>	paper rings with allergen	01 x 10
EP WD25-0110	22861	<b>WD25: Red Cedar</b>	paper rings with allergen	01 x 10
EP WS01-0125	22167	<b>WS1: Weeds I (W1, W6, W9, W10, W11)</b>	paper rings with allergen mixture	01 x 25
EP WS13-0125	22168	<b>WS13: Weeds II (W6, W9, W10, W12)</b>	paper rings with allergen mixture	01 x 25
EP WS14-0110	22965	<b>WS14: Flowers I (W30, W25, W33, W40)</b>	paper rings with allergen mixture	01 x 10
EP WS15-0110	22973	<b>WS15: Flowers II (W51, W22, W23, W7)</b>	paper rings with allergen mixture	01 x 10
EP WS16-0110	22695	<b>WS16: Weeds III (W1, W6, W7, W8, W12)</b>	paper rings with allergen mixture	01 x 10
EP WS17-0110	22699	<b>WS17: Weeds IV (W9, W10, W11, W100)</b>	paper rings with allergen mixture	01 x 10